

Product datasheet for RC226393L3V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

MET (NM 001127500) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: MET (NM_001127500) Human Tagged ORF Clone Lentiviral Particle

Symbol:

AUTS9; c-Met; DFNB97; HGFR; RCCP2 Synonyms:

Mammalian Cell

Selection:

Puromycin

pLenti-C-Myc-DDK-P2A-Puro (PS100092) Vector:

Myc-DDK Tag:

NM 001127500 ACCN:

ORF Size: 4224 bp

ORF Nucleotide

Sequence:

The ORF insert of this clone is exactly the same as(RC226393).

The molecular sequence of this clone aligns with the gene accession number as a point of OTI Disclaimer: reference only. However, individual transcript sequences of the same gene can differ through

naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: NM 001127500.1

RefSeq ORF: 4227 bp Locus ID: 4233

UniProt ID: P08581 **Cytogenetics:** 7q31.2

Protein Families: Druggable Genome, Protein Kinase, Transmembrane

Protein Pathways: Adherens junction, Axon guidance, Colorectal cancer, Cytokine-cytokine receptor interaction,

Endocytosis, Epithelial cell signaling in Helicobacter pylori infection, Focal adhesion,

Melanoma, Pathways in cancer, Renal cell carcinoma







MW:

157.71 kDa

Gene Summary:

This gene encodes a member of the receptor tyrosine kinase family of proteins and the product of the proto-oncogene MET. The encoded preproprotein is proteolytically processed to generate alpha and beta subunits that are linked via disulfide bonds to form the mature receptor. Further processing of the beta subunit results in the formation of the M10 peptide, which has been shown to reduce lung fibrosis. Binding of its ligand, hepatocyte growth factor, induces dimerization and activation of the receptor, which plays a role in cellular survival, embryogenesis, and cellular migration and invasion. Mutations in this gene are associated with papillary renal cell carcinoma, hepatocellular carcinoma, and various head and neck cancers. Amplification and overexpression of this gene are also associated with multiple human cancers. [provided by RefSeq, May 2016]