

## Product datasheet for **RC226417L4V**

### **BRG1 (SMARCA4) (NM\_001128844) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	BRG1 (SMARCA4) (NM_001128844) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SMARCA4
Synonyms:	BAF190; BAF190A; BRG1; CSS4; hSNF2b; MRD16; RTPS2; SNF2; SNF2-beta; SNF2L4; SNF2LB; SWI2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001128844
ORF Size:	4942 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226417).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001128844.1</a> , <a href="#">NP_001122316.1</a>
RefSeq Size:	5779 bp
RefSeq ORF:	4944 bp
Locus ID:	6597
UniProt ID:	<a href="#">P51532</a>
Cytogenetics:	19p13.2
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

**MW:** 185.1 kDa

**Gene Summary:** The protein encoded by this gene is a member of the SWI/SNF family of proteins and is similar to the brahma protein of *Drosophila*. Members of this family have helicase and ATPase activities and are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein is part of the large ATP-dependent chromatin remodeling complex SNF/SWI, which is required for transcriptional activation of genes normally repressed by chromatin. In addition, this protein can bind BRCA1, as well as regulate the expression of the tumorigenic protein CD44. Mutations in this gene cause rhabdoid tumor predisposition syndrome type 2. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2012]