

Product datasheet for **RC226444**

NACA (NM_001113203) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NACA (NM_001113203) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NACA
Synonyms:	HSD48; NAC-alpha; NACA1; skNAC
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC226444 representing NM_001113203
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCCCGCGAAGCCACAGAAACCGTCCCTGCTACAGAGCAGGAGTTGCCGCAGCCCCAGGCTGAGACAG
CTGTGCTACCTATGTCTTCAGCCTTGAGTGTCACTGCTGCCTTAGGGCAGCCTGGACCTACCCTCCCCC
TCCTTGCTCTCCTGCCCAACAGTGCCTCTCTCAGCTGCTAACCCAGGCTTCCCCATTCCCTTCCCC
TCTACTATTGCCTCGACCCCTTTAGAAGTTCCTTTTCCCAGTCATCCTCTGGAACAGCCCTACCTTTGG
GAACTGCCCTGAAGCCCCAACCTTCTACCAACCTAATAGGGCCTCCCATCTCCCAGCTGCCTTAGC
TCTAGCCTCTCCCATGATAGCTCCAACCTGAAAGGGACCCCTTCTCTCAGCTCCCTAGCTCTGGTT
GCCCTGGTCCCCACTCAGTTCAGAAGAGTTCTGCTTTTCCACCTAACCTTCTTACTTCACTCCTTCAG
TGGCTGTAGCTGAGTCAGGTCAGTGATAACTCTGTGAGCTCCCATGCTCCCTCAGAACCAAAGACTAA
TCTTAATAAAGTTCCTCTGAGGTAGTCCCTAATCCAAAAGGCACCCCAAGCCCTCCATGTATAGTCAGT
ACTGTTCCCTTACCACTGTGTGACTCCCATGGCCTCTATCAATCTGGAGTGGCCTCCCTTCTCAGACAA
CACCCACAACCTACCTAGCCATCGCTTCCCTCAAGTCAAAGATACCACCATTTCTCAGTTCGTATTTTC
TCCAAAAACCCAGGAAGCCTCAGCCTGAAGGGCCTGTAGTCCACCTGCTGCCTTATCTTTTCAACT
CAGTCTCTTCTGTGGTACCTCTTCTCAAAGACTGCGGGTCCCAACACCCCCAGATTTTCCATTT
CTCTGGGCTCTCATCTTGACCTTTACATCAGAGTTCTTTTGGTCTGTCCAACCTTTAGGTCAAACAGG
TCCTAGTGCTTTGTCAGACCCTACAGTGAAGACCATTTCTGTAGATCATTCTCCACAGGGGCTCTTAT
CCTTCTCAGAGATCTGTAATTCCTCCCTTCTTCCAGAAATGAGGTAGTCTGCTACTGCAGCCCTCTGG
TTCCAGTGGTGGCTCCATCTGTTGACAAAGGTCCTCTACCATCTCTAGCATAACCTGCAGCCCTCTGG
CTCCTTAAATGTAGCTACCTCTTTTTCATTATCTCCTACAACTCTCTCATTCTCAAAGCTCTCCTAAT
GCCACTTATCATTATCCTTTAGTGGCCCAATGCCGTTTCTTCTGTTGGAACCCCACTTGTGGTGA
CTAACCCCTGTACAATTGCTGCAGCACCTACTACTACCTTTGAGGTAGCTACTTGTGTTTCTCCTCAAT
GTCATCAGGTCCCATAGTAACATAGAACCAACTTCCCTGCTGCCTTGGTTATGGCACCTGTGGCTCCC
AAAGAGCCTTCTACTCAAGTAGCAACCACTCTGAGGATACCAGTCTCTCCTCTGCCAGACCCTGAAG
ACCTCAAAAATCTCCCCAGTTCAGTATTGGTTAAATTTCCAACACAAAAAGACCTCAAACCTGTACCTGC
CTCTCTGAAGGAGCCCTTTCTCTCCAGCCCAAGCAGGACTCACCAAGAAAGACCCTACTGTATTA
CGTTAGTCCAGGAGCCCTAAAATTCCTTCTTTCAAAGTACATCCTCTTCTCCAGAGATACCTC
TTTCTCCTGAAGCCACCCTAGCAAAGAAAAGCCTTGGGGAGCCTCTCCCTATAGTGGCTGCATTTCCCTT
GAAAAGTGTGACCTGCCGGGGTGGCTCCCAACTGCCAAAGCAGCTGCCTTTGAGAAGTCTTCTCT
AAACCTGAATCAGCATCTGTCTCTGCAGCACCTCCCCACCAGTCTCTGCTCCTTGTCTCCCTCCCCAG
TTCCCACTCTGCCTCCTAAACAGCAATTTCTGCCGTCCTCTCCTGGGCTGGTGTGGAATCACCTCTAA
ACCCCTTGCCCTGTGATGAGGATGAGCTGCTGCCTCTGATTCCCCCGAACCAATCTCTGGGGGAGTG
CCTTCCAGTCGGTCTCGTCAACATGCCACCCCTAAATCTGCTGGAATCCCTGTCCCAACCCCTCTG
CCAAGCAACCTGTTACGAAGAACAACAAGGGTCTGGAACAGAATCTGACAGTGAATCAGTACCAGA
GCTTGAAGAACAGGATTCACCCAGGCAACCACACAACAAGCCAGCTGGCGGCAGCAGCTGAAATTGAT
GAAGAACCAGTCAGTAAAGCAAAACAGAGTCGGAGTGAAAAGAAGGCACGGAAGGCTATGTCCAACTGG
GTCTTCGCGAGGTTACAGGAGTTACTAGAGTCACTATCCGAAATCTAAGAATATCCTCTTTGTATCAC
AAAACCAGATGTCTACAAGAGCCCTGCTTCAAGTACTTACATAGTTTTTGGGAAGCCAAGATCGAAGAT
TTATCCAGCAAGCACAACCTAGCAGCTGCTGAGAAATTCAAAGTCAAGGTGAAGCTGTCTCAAACATTC
AAGAAAACACACAGACTCCAACCTGTACAAGAGGAGAGTGAAGAGGAAGAGGTCGATGAAACAGGTGAGA
AGTTAAGGACATTGAATTGGTCAATGTACAAGCAAATGTGTCGAGAGCAAAGGCAGTCCGAGCCCTGAAG
AACACAGTAATGATATTGTAATGCGATTATGGAATTAACAATG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC226444 representing NM_001113203
 Red=Cloning site Green=Tags(s)

MPGEATETVPATEQELPQPQAETA V L P M S S A L S V T A A L G Q P G P T L P P P C S P A P Q Q C P L S A A N Q A S P F P S P
 S T I A S T P L E V P F P Q S S S G T A L P L G T A P E A P T F L P N L I G P P I S P A A L A L A S P M I A P T L K G T P S S S A P L A L V
 A L A P H S V Q K S S A F P P N L L T S P P S V A V A E S G S V I T L S A P I A P S E P K T N L N K V P S E V V P N K G T P S P P C I V S
 T V P Y H C V T P M A S I Q S G V A S L P Q T T P T T L A I A S P Q V K D T T I S S V L I S P Q N P G S L S L K G P V S P P A A L S L S T
 Q S L P V V T S S Q K T A G P N T P P D F P I S L G S H L A P L H Q S S F G S V Q L L G Q T G P S A L S D P T V K T I S V D H S S T G A S Y
 P S Q R S V I P P L P S R N E V P A T V A A F P V V A P S V D K G P S T I S S I T C S P S G S L N V A T S F S L S P T T S L I L K S S P N
 A T Y H Y P L V A Q M P V S S V G T T P L V V T N P C T I A A A P T T T F E V A T C V S P P M S S G P I S N I E P T S P A A L V M A P V A P
 K E P S T Q V A T T L R I P V S P P L P D P E D L K N L P S S V L V K F P T Q K D L Q T V P A S L E G A P F S P A Q A G L T T K D P T V L
 P L V Q A A P K N S P S F Q S T S S S P E I P L S P E A T L A K K S L G E P L P I V A A F P L E S A D P A G V A P T T A K A A A F E K V L P
 K P E S A S V S A A P S P P V S L P L A P S P V P T L P P K Q Q F L P S S P G L V L E S P S K P L A P A D E D E L L P L I P P E I S G G V
 P F Q S V L V N M P T P K S A G I P V P T P S A K Q P V T K N N K G S G T E S D S D E S V P E L E E Q D S T Q A T T Q Q A Q L A A A A E I D
 E E P V S K A K Q S R S E K K A R K A M S K L G L R Q V T G V T R V T I R K S K N I L F V I T K P D V Y K S P A S D T Y I V F G E A K I E D
 L S Q Q A Q L A A A E K F K V Q G E A V S N I Q E N T Q T P T V Q E E S E E E E V D E T G V E V K D I E L V M S Q A N V S R A K A V R A L K
 N N S N D I V N A I M E L T M

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:




ACCN: NM_001113203

ORF Size: 2775 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001113203.3
RefSeq Size:	3191 bp
RefSeq ORF:	2778 bp
Locus ID:	4666
UniProt ID:	Q13765
Cytogenetics:	12q13.3
Protein Families:	Druggable Genome, Transcription Factors
MW:	94.7 kDa
Gene Summary:	This gene encodes a protein that associates with basic transcription factor 3 (BTF3) to form the nascent polypeptide-associated complex (NAC). This complex binds to nascent proteins that lack a signal peptide motif as they emerge from the ribosome, blocking interaction with the signal recognition particle (SRP) and preventing mistranslocation to the endoplasmic reticulum. This protein is an IgE autoantigen in atopic dermatitis patients. Alternative splicing results in multiple transcript variants, but the full length nature of some of these variants, including those encoding very large proteins, has not been determined. There are multiple pseudogenes of this gene on different chromosomes. [provided by RefSeq, Feb 2016]

Product images:



Circular map for RC226444