

Product datasheet for RC226508

RAB34 (NM 001144943) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: RAB34 (NM_001144943) Human Tagged ORF Clone

Tag: Myc-DDK
Symbol: RAB34

Synonyms: NARR; RAB39; RAH

Mammalian Cell Neomycin

Selection:

Vector:pCMV6-Entry (PS100001)E. coli Selection:Kanamycin (25 ug/mL)

ORF Nucleotide >RC226508 representing NM_001144943

Sequence: Red=Cloning site Blue=ORF Green=Tags(s)

 ${\tt TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC}$

GCCGCGATCGCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com Protein Sequence: >RC226508 representing NM_001144943

Red=Cloning site Green=Tags(s)

MSHLPGLELRREAPPLLGPLLSPFPLPAGSWHRQMLRSSLRFPITNSAGAPCKAAGRMNILAPVRRDRVL AELPQCLRKEAALHGHKDFHPRVTCACQEHRTGTVGFKISKVIVVGDLSVGKTCLINRFCKDTFDKNYKA TIGVDFEMERFEVLGIPFSLQLWDTAGQERFKCIASTYYRGAQAIIIVFNLNDVASLEHTKQWLADALKE NDPSSVLLFLVGSKKDLSTPAQYALMEKDALQVAQEMKAEYWAVSSLTGENVREFFFRVAALTFEANVLA

ELEKSGARRIGDVVRINSDDSNLYLTASKKKPTCCP


TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8051 d08.zip

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM 001144943

ORF Size: 948 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:

- 1. Centrifuge at 5,000xg for 5min.
- 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: <u>NM 001144943.1</u>, <u>NP 001138415.1</u>

 RefSeq ORF:
 951 bp

 Locus ID:
 83871

 Cytogenetics:
 17q11.2

Protein Families: Druggable Genome

MW: 35 kDa

Gene Summary: This gene encodes a protein belonging to the RAB family of proteins, which are small GTPases


involved in protein transport. This family member is a Golgi-bound member of the secretory

pathway that is involved in the repositioning of lysosomes and the activation of

macropinocytosis. Alternative splicing of this gene results in multiple transcript variants. An alternatively spliced transcript variant produces the nine-amino acid residue-repeats (NARR) protein, which is a functionally distinct nucleolar protein resulting from a different reading


frame. [provided by RefSeq, Dec 2016]


Product images:



Circular map for RC226508







Western blot validation of overexpression lysate (Cat# [LY428595]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC226508 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).

Coomassie blue staining of purified RAB34 protein (Cat# [TP326508]). The protein was produced from HEK293T cells transfected with RAB34 cDNA clone (Cat# RC226508) using MegaTran 2.0 (Cat# [TT210002]).