

## Product datasheet for **RC226526L3V**

### **C17orf104 (MEIOC) (NM\_001145080) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	C17orf104 (MEIOC) (NM_001145080) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C17orf104
Synonyms:	C17orf104
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001145080
ORF Size:	2856 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226526).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001145080.1</a>
RefSeq ORF:	2859 bp
Locus ID:	284071
UniProt ID:	<a href="#">A2RUB1</a>
Cytogenetics:	17q21.31
MW:	107.4 kDa



[View online »](#)

**Gene Summary:**

Is required for meiosis completion in both male and female germ cells. Confers stability to numerous meiotic mRNAs in gonads allowing proper initiation and progression into meiosis prophase I. The function may involve YTHDC2 and is independent of induction by retinoic acid (RA). Maintains an extended meiotic prophase I by properly promoting the transition from a mitotic to a meiotic cell cycle program by binding transcripts through its interaction with YTHDC2 that regulate the mitotic cell cycle.[UniProtKB/Swiss-Prot Function]