

Product datasheet for **RC226571A1V**

Human IL11RA (NM_001142784) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human IL11RA (NM_001142784) AAV Particle
Tag:	Myc-DDK
Symbol:	IL11RA
Synonyms:	CRSDA
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC226571 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGCAGCAGCTGCTCAGGGCTGAGCAGGGTCTGGTGGCCGTGGCTACAGCCCTGGTGTCTGCCTCCT
CCCCCTGCCCCAGGCCCTGGGGCCCCCAGGGGTCCAGTATGGGCAGCCAGGCAGGTCCGTGAAGCTGTG
TTGTCTGGAGTGACTGCCGGGACCCAGTGTCTGGTTTCGGGATGGGGAGCCAAAGCTGCTCCAGGGA
CCTGACTCTGGGCTAGGGCATGAACTGGTCTGGCCAGGCAGACAGCACTGATGAGGGACCTACATCT
GCCAGCCCTGGATGGTGCACCTGGGGGCACAGTGACCCTGCAGCTGGGCTACCCTCCAGCCCGCCCTGT
TGTCTCTGCCAAGCAGCCGACTATGAGAACTTCTCTTGCACTTGGAGTCCAGCCAGATCAGCGGTTTA
CCCACCCGCTACCTCACCTCCTACAGGAAGAAGACAGTCTAGGAGCTGATAGCCAGAGGAGGAGTCCAT
CCACAGGGCCCTGGCCATGCCCACAGGATCCCCTAGGGGCTGCCCGTGTGTTGTCCACGGGGCTGAGTT
CTGGAGCCAGTACCGGATTAATGTGACTGAGGTGAACCCACTGGGTGCCAGCACAGCCTGCTGGATGTG
AGCTTGCAGAGCATCTTGCGCCCTGACCCACCCAGGGCTGCCGGTAGAGTCAGTACCAGGTTACCCCC
GACGCTGCGAGCCAGCTGGACATACCCTGCCTCCTGGCCGTGCCAGCCCACTTCTGCTCAAGTTCCG
TTTGAGTACCGTCCGGCGCAGCATCCAGCCTGGTCCACGGTGGAGCCAGCTGGACTGGAGGAGGTGATC
ACAGATGCTGTGGCTGGGCTGCCCCATGCTGTACGAGTCAGTGCCCGGACTTTCTAGATGCTGGCACCT
GGAGCACCTGGAGCCCGAGGCCTGGGGAATCCGAGCACTGGGACCATACCAAAGGAGATACCAGCATG
GGGCCAGCTACACACGACCCAGAGGTGGAGCCTCAGGTGGACAGCCCTGCTCCTCCAAGGCCCTCCCTC
CAACCACACCCTCGGCTACTTGATCACAGGGACTCTGTGGAGCAGGTAGCTGTGCTGGCGTCTTTGGGAA
TCCTTTCTTCTGGGACTGGTGGCTGGGGCCCTGGCACTGGGGCTCTGGCTGAGGCTGAGACGGGGTGG
GAAGGATGGATCCCCAAGCCTGGGTTCTTGGCCTCAGTGATTCCAGTGGACAGGCGTCCAGGAGCTCCA
AACCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC226571 protein sequence Red=Cloning site Green=Tags(s) MSSSCSGLSRVLVAVATALVSASSPCPQAWGPPGVQYGQPGRSVKLCCPGVTAGDPVSWFRDGEPKLLQG PDSGLGHEL VLAQADSTDEGTICQTLDGALGGTVTLQLGYPPARPVVSCQAADYENFCTWSPSQISGL PTRYLT SYRKKT VLGADSQRRSPSTGPWPCPDPLGAARCVVHGAEFWSQYRINVTEVNPLGASTRLLDV SLQSILRPDPPQGLRVESVPGYPRRLRASWTYPASWPCQPHFLKFRLLQYRPAQHPAWSTVEPAGLEEV TDAVAGLPHAVRVSARDFLDAGTWSTWSPEAWGTPSTGTIPKEIPAWGQLHTQPEVEPQVDSAPPRPSL QPHPRLLDHRDSVEQVAVLASLGILSFLGLVAGALALGLWRLRRGGKDGSPKPGFLASVIPVDRRPGAP NL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001142784
ORF Size:	1266 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001142784.1
RefSeq Size:	1728 bp
RefSeq ORF:	1269 bp
Locus ID:	3590
UniProt ID:	Q14626
Cytogenetics:	9p13.3
MW:	45.2 kDa