

Product datasheet for **RC226577L4V**

SAR1 (SAR1A) (NM_001142648) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | SAR1 (SAR1A) (NM_001142648) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | SAR1 |
| Synonyms: | masra2; SAR1; Sara; SARA1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001142648 |
| ORF Size: | 594 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC226577). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001142648.1 |
| RefSeq Size: | 3088 bp |
| RefSeq ORF: | 597 bp |
| Locus ID: | 56681 |
| UniProt ID: | Q9NR31 |
| Cytogenetics: | 10q22.1 |
| MW: | 22.4 kDa |



[View online »](#)

Gene Summary:

Involved in transport from the endoplasmic reticulum to the Golgi apparatus (By similarity). Required to maintain SEC16A localization at discrete locations on the ER membrane perhaps by preventing its dissociation. SAR1A-GTP-dependent assembly of SEC16A on the ER membrane forms an organized scaffold defining endoplasmic reticulum exit sites (ERES). [UniProtKB/Swiss-Prot Function]