

Product datasheet for **RC226579A1V**

Human NONO (NM_001145409) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human NONO (NM_001145409) AAV Particle
Tag:	Myc-DDK
Symbol:	NONO
Synonyms:	MRXS34; NMT55; NRB54; P54; P54NRB; PPP1R114
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC226579 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCAGAGTAATAAACTTTTAACTTGGAGAAGCAAAACCATACTCCAAGAAAGCATCATCAACATCACC
 ACCAGCAGCAGCACCACCAGCAGCAACAGCAGCAGCCGCCACCACCGCCAATACCTGCAAATGGGCAACA
 GGCCAGCAGCCAAAATGAAGGCTTGACTATTGACCTGAAGAATTTAGAAAACAGGAGAGAAGACCTTC
 ACCCAACGAAGCCGTCCTTTTGTGGAAATCTTCTCCGACATCACTGAGGAAGAAATGAGGAACTAT
 TTGAGAAATATGGAAAGGCAGGCGAAGTCTTCATTATAAGGATAAAGGATTTGGCTTTATCCGCTTGA
 AACCCGAACCCTAGCGGAGATTGCCAAGTGGAGCTGGACAATATGCCACTCCGTGGAAGCAGCTGCGT
 GTGCGCTTTGCTGCCATAGTCATCCCTACAGTTCGAAACCTTCTCAGTATGTGTCAACGAAGTGC
 TGAAGAAGCCTTTTCTGTGTTGGCCAGGTAGAGAGGCTGTAGTCATTGTGGATGATCGAGGAAGGCC
 CTCAGGAAAAGCATTGTTGAGTCTCAGGGAAGCCAGCTGCTCGAAAGCTCTGGACAGATGCAGTGAA
 GGCTCCTTCTGCTAACCACATTTCTCTGTCCTGTGACTGTGGAGCCATGGACCAGTTAGATGATGAAG
 AGGGACTTCCAGAGAAGCTGGTTATAAAAAACCAGCAATTTACAAGGAACGAGAGCAGCCACCCAGATT
 TGCACAGCCTGGCTCCTTTGAGTATGAATATGCCATGCGCTGGAAGGCACTCATTGAGATGGAGAAGCAG
 CAGCAGGACCAAGTGGACCGCAACATCAAGGAGGCTCGTGAGAAGCTGGAGATGGAGATGGAAGCTGCAC
 GCCATGAGCACCAGGTGCTAATGAGACAGGATTTGATGAGGCGCCAAGAAGAACTTCGGAGGATGGA
 AGAGCTGCACAACCAAGAGGTGCAAAAACGAAAGCAACTGGAGCTCAGGAGGAGGAAGAGCGCAGGCGC
 CGTGAAGAAGAGATGCGGCGGCAGCAAGAAGAAATGATGCGGCGACAGCAGGAAGGATTCAGGGAACTT
 TCCCTGATGCGAGAGAGCAGGAGATTCGGATGGGTGAGTGGCTATGGGAGGTGCTATGGGCATAACAA
 CAGAGGTGCCATGCCCTGCTCTGTGCCAGCTGGTACCCAGCTCCTCCAGGACCTGCCACTATGATG
 CCGGATGGAACCTTTGGGATTGACCCACCAACAACCTGAACGCTTTGGTCAGGCTGCTACAATGGAAGGAA
 TTGGGGCAATTGGTGGAACCTCTCTGCATTCAACCGTCAGCTCCTGGAGCTGAATTTGCCCAACAA
 ACGTCGCCGATAC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC226579 protein sequence
 Red=Cloning site Green=Tags(s)

MQSNKTFNLEKQNHTRPKHHQHQQHQQHQQHQQPPPPPIPANGQQASSQNEGLTIDLKNFRKPGEKTF
 TQRSRLFVGNLPPDITEEEMRKLFEKYKGAGEVFIHKDKGFGFIRLETRTLAEIAKVELDNMPLRGKQLR
 VRFACHSASLTVRNLPQYVSNELLEEAFSVFGQVERAVVIVDDRGRPSGKGIVEFSGKPAARKALDRCSE
 GSFLLTFFRPVTVPEMDQLDDEEGLPEKLVIKNQQFHKEREQPPRFAQPGSFYEYAMRWKAL IEMEKQ
 QQDQVDRNIKEAREKLEMEMEAAARHEHQVLMRQDLMRQEELRRMEELHNQEVQKRQLELRQEEERRR
 REEEMRRQQEEMRRQEGFKGTFPDAREQEIRMGQMAMGGAMGINNRGAMPAPVPAGTPAPPGPATMM
 PDGTLGLTPPTTERFQGAATMEGIGAIGTPPAFNRAAPGAEPNKRRTY

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species:

Human

Serotype:

AAV-2

ACCN:

NM_001145409

ORF Size:

1413 bp

Buffer:

PBS with 0.001% Pluronic F68

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001145409.1
RefSeq Size:	3057 bp
RefSeq ORF:	1416 bp
Locus ID:	4841
UniProt ID:	Q15233
Cytogenetics:	Xq13.1
MW:	54.2 kDa