

## Product datasheet for **RC226664**

### EIF5A (NM\_001143762) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** EIF5A (NM\_001143762) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** EIF5A  
**Synonyms:** eIF-4D; EIF-5A; EIF5A1; eIF5AI  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RC226664 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCAGATGACTTGGACTTCGAGACAGGAGATGCAGGGGCTCAGCCACCTTCCAATGCAGTGCTCAG  
CATTACGTAAGAATGGCTTTGTGGTGCTCAAAGGCCGGCCATGTAAGATCGTCGAGATGTCTACTTCGAA  
GACTGGCAAGCACGGCCACGCCAAGGTCCATCTGGTTGGTATTGACATCTTTACTGGGAAGAAATATGAA  
GATATCTGCCGTCAACTATAATATGGATGTCCCAACATCAAAGGAATGACTCCAGCTGATTGGCA  
TCCAGGATGGGTACCTATCACTGCTCCAGGACAGCGGGGAGGTACGAGAGGACCTTCGTCCTCCCTGAGGG  
AGACCTTGCCAAGGAGATTGAGCAGAAGTACGACTGTGGAGAAGAGATCCTGATCACGGTGTCTGCC  
ATGACAGAGGAGGCAGCTGTTGCAATCAAGGCCATGGCAAAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC226664 protein sequence  
Red=Cloning site Green=Tags(s)

MADDLDFETGDAGASATFPMQCSALRKNFVVLKGRPCKIVEMSTSKTGKHGKAVHLVGDIDFTGKKYE  
DICPSTHNMDVPIKRNDQLIGIQDGYLSLLQDSGEVREDLRLPEGDLGKEIEQKYDCGEEILITVLSA  
MTEEAAVAIAKAMAK

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk6053\\_g11.zip](https://cdn.origene.com/chromatograms/mk6053_g11.zip)



[View online >](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM\_001143762

ORF Size: 462 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001143762.1](#), [NP\\_001137234.1](#)

RefSeq Size: 1340 bp


RefSeq ORF: 465 bp

Locus ID: 1984


UniProt ID: [P63241](#)  
 Cytogenetics: 17p13.1  
 MW: 16.8 kDa

**Gene Summary:** mRNA-binding protein involved in translation elongation. Has an important function at the level of mRNA turnover, probably acting downstream of decapping. Involved in actin dynamics and cell cycle progression, mRNA decay and probably in a pathway involved in stress response and maintenance of cell wall integrity. With syntenin SDCBP, functions as a regulator of p53/TP53 and p53/TP53-dependent apoptosis. Regulates also TNF-alpha-mediated apoptosis. Mediates effects of polyamines on neuronal process extension and survival. May play an important role in brain development and function, and in skeletal muscle stem cell differentiation. Also described as a cellular cofactor of human T-cell leukemia virus type I (HTLV-1) Rex protein and of human immunodeficiency virus type 1 (HIV-1) Rev protein, essential for mRNA export of retroviral transcripts.[UniProtKB/Swiss-Prot Function]


**Product images:**



Circular map for RC226664



Western blot validation of overexpression lysate (Cat# [LY428327]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC226664 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified EIF5A protein (Cat# [TP326664]). The protein was produced from HEK293T cells transfected with EIF5A cDNA clone (Cat# RC226664) using MegaTran 2.0 (Cat# [TT210002]).