

Product datasheet for **RC226682L3V**

DAGLB (NM_001142936) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DAGLB (NM_001142936) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DAGLB
Synonyms:	DAGLBETA; KCCR13L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001142936
ORF Size:	1629 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226682).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001142936.1 , NP_001136408.1
RefSeq ORF:	1632 bp
Locus ID:	221955
UniProt ID:	Q8NCG7
Cytogenetics:	7p22.1
Protein Families:	Transmembrane
MW:	59.7 kDa



[View online »](#)

Gene Summary:

Catalyzes the hydrolysis of diacylglycerol (DAG) to 2-arachidonoyl-glycerol (2-AG), the most abundant endocannabinoid in tissues. Required for axonal growth during development and for retrograde synaptic signaling at mature synapses.[UniProtKB/Swiss-Prot Function]