

Product datasheet for **RC226690**

FAM111B (NM_001142704) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM111B (NM_001142704) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	FAM111B
Synonyms:	CANP; POIKTMP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC226690 representing NM_001142704
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAATCCATGAAGACTGAAGAAAACAAGTCATTTAGCGCTATGGAAGATGACCAGAGGACTAGACCTG
 AAGTTTCAAAGGATACTGTCATGAAGCAGACACATGCTGACACACCTGTTGATCATTGTCTATCTGGCAT
 AAGAAAGTGTAGCAGCACCTTTAAGCTTAAAGTGAAGTCAACAAGCATGAAACAGCCCTTGAAATGCAG
 AATCCAAATTTGAACAATAAAGAATGTTGTTTACCTTTACGTTGAATGGAACTCCAGAAAATTAGACC
 GTAGTGTGTTTACAGCATATGGTAAACCCAGCGAGAGTATCTACTCAGCCCTGAGTGCTAATGACTATTT
 CAGTGAAGGATAAAGAATCAGTTTAATAAGAACATTATTGTTTATGAAGAAAAGACAATAGATGGACAT
 ATAAATTTAGGAATGCCTCTCAAGTGCCTGCTAGTGATTCTCATTTTAAATTTACATTTGGTCAAAGAA
 AGAGTAGCAAAGAAGATGGACACATATTACGCCAATGTGAAAATCCAAACATGGAATGCATTCTTTTTCA
 TGTTGTTGCTATAGGAAGGACAAGAAAGAAGATTGTTAAGATCAACGAACCTCATGAAAAAGGAAGTAA
 CTTTGTATTTATGCCTTGAAAGGGTGAGACTATTGAAGGAGCCTTATGCAAGGATGGCCGTTTTCGGTCTG
 ACATAGGTGAATTTGAATGGAACTAAAGGAAGGTCATAAGAAAAATTTATGGAAAACAGTCCATGGTGGA
 TGAAGTATCTGGAAAAGTCTTAGAAATGGACATTTCAAAAAAAGCATTACAACAGAAAGATATCCAT
 AAAAAATTAACAGAATGAAAGTGCCACTGATGAAATTAATCACCAGAGTCTGATACAGCTAAGAAAA
 AAGTCCACAAACCAAGAAAGATGGAGAGACCAAGATGTAGAACACAGCAGAGAGCAAATTCCTCCACC
 TCAGGATCTAAGCCATTATTTAAAGATAAACTCGCCAGACAATTCCAGGATTAGAAATTTACTTTT
 TGTAGTTTGCCCGAAAAATATAGGCAATAAACTCACAGTTAGACGGAGGCCGCATCTGGGTAGGCGGT
 ATGCTATTAATCTGGATGTCCAAAGGAGGCAATTAATCTCTTAAGAATTATCAAACGTGAATGAAGC
 CATAATGCATCAGTATCCGAATTTTAAAGAGGAGGCACAGTGGGTAAAGAAAATATTTTCGGGAAGAACA
 AAGAGAATGAATCTTTCACCAGCTAAGCAATTCACATATATAAAAGGACTTCGAAAAATGACTGCAA
 ATTCTGTTTCAGTTGCAACCTGCGAACAGCTTACATATTATAGCAAGTCAGTTGGGTTTCATGCAATGGGA
 CAATAATGGAAACACAGGTAATGCTACTTGCTTTGTCTTCAATGGTGGTTATATTTTACCTGTCGACAT
 GTTGTACATCTTATGGTGGGTAAAAACACACATCCAAGTTTGTGGCCAGATATAATTAGCAAATGTGCGA
 AGGTAACCTTCACTTATACAGAGTTCTGCCCTACTCCTGACAATTGGTTTTCCATTGAGCCATGGCTTAA
 AGTGTCCAATGAAATCTAGATTATGCCATTTTAAACTAAAAGAAAATGGAATGCGTTTCCTCCAGGA
 CTATGGCGACAGATTTCTCCTCAACCATCTACTGGTTTGATTTATTTAATTGGTCATCCTGAAGGCCAGA
 TCAAGAAAATAGATGGTTGTACTGTGATTCTCTAAACGAACGATTGAAAAAATATCCAAACGATTGTCA
 AGATGGGTTGGTAGATCTCTATGATACCACAGTAATGTATACTGTATGTTTACCCAAAGAAGTTTCCTA
 TCAGAGGTTTGGAAACACACACAGCTTAGTTATGATACTTGTTTCTCTGATGGGTCTCAGGCTCCCCAG
 TGTTTAAATGCATCTGGCAAAATGGTTGCTTTCATACCTTTGGGCTTTTTATCAACGAGGATTTAATGT
 GCATGCCCTTATTGAATTTGGTTATTCTATGGATTCTATTCTTTGTGATATTAAGACAAATGAGAGC
 TTGTATAAATCATTAAATGATGAGAACTTGAGACCTACGATGAAGAGAAAGGTAACAAGAGTCATCAC
 TTCAAGATCATCAGATTGAACCCATGGAATGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC226690 representing NM_001142704
Red=Cloning site Green=Tags(s)

MNSMKT EENKSF SAMEDDQRT RPEVSKD TVMKQTHADTPVDHCLSGIRKCSSTFKLKSEVNKHETALEMQ
NPNLNNKECCFTFTLNGNSRKLD RSVFTAYGKPSEIYSAL SANDYF SERIKNQFNKNIIVYEKTI DGH
INLGMPLKCLPSDSHFKITFGQRKSSKEDGHILRQCENPNMECILFHVVAIGRTRKKIVKINELHEKGSK
LCIYALKGETIEGALCKDGRFRSDIGEFEWLKEGHHKIYGQKSMVDEVSGKVL EMDISKKALQQKDIH
KKIKQNESATDEINHQSLIQSKKKVHKPKKDGETKDVHSREQILPPQDLSHYIKDKTRQTIPIRIRNYF
CSLPRKYRQINSQVRRRPHLGRRYAINL DVQKEAINLLKNYQTLNEAIMHQYPNFKEEAQWVRKYFREEQ
KRMNLSPAKQFNIYKDFGKMTANSVSATCEQLTYYSKSVGFMQWDNNGNTGNATCFVNGGYIFTCRH
VVHLMVGKNTHPSLWPDII SKCAKVTFYTFECPTPDNWF SIEPWLKVSNNENLDYAILKLKENGNAFFPG
LWRQISQPSTGLIYLIGHPEGQIKIDGCTVALPTNERLKKYPNDCQDGLVDLYDTTSNVYCMFTQRSFL
SEVWNHTLSYDTCFSDGSSQSPVFNASGKLVALTFLGYQRGFNVHALIEFGYSMDSILCDIKKTNES
LYKSLNDEKLEITYDEEKGQOESSLODHOT EPMEC

TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM 001142704

ORF Size: 2202 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001142704.1](#), [NP_001136176.1](#)

RefSeq Size: 3362 bp

RefSeq ORF: 2115 bp

Locus ID: 374393


UniProt ID: [Q6SJ93](#)

Cytogenetics: 11q12.1

MW: 84.7 kDa

Gene Summary: This gene encodes a protein with a trypsin-like cysteine/serine peptidase domain in the C-terminus. Mutations in this gene are associated with an autosomal dominant form of hereditary fibrosing poikiloderma (HFP). Affected individuals display mottled pigmentation, telangiectasia, epidermal atrophy, tendon contractures, and progressive pulmonary fibrosis. Alternative splicing results in multiple transcript variants encoding distinct isoforms. A paralog of this gene which also has a trypsin-like peptidase domain, FAM111A, is located only 16 kb from this gene on human chromosome 11q12.1. [provided by RefSeq, Apr 2014]

Product images:



Circular map for RC226690