

Product datasheet for **RC226699**

HNRPLL (HNRNPLL) (NM_001142650) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	HNRPLL
Synonyms:	HNRPLL; SRRF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >RC226699 representing NM_001142650
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGTTGGCTCGGGAGACGTACGAGGAGGACCGGGAGTACGAGAGCCAGGCCAAGCGTCTCAAGACCGAGG
 AGGGGGAGATCGACTACTCGGCCGAGGAAGGCGAGAACC CGCGGAAGCGACGCCCCGGGGCGGGGCGA
 TGGCGGCGGGCGGCCGGAGCTTCTCTCAGCCGGAGGCGAGGTGGAAGTCATCAAAAGTTTCTGTTTCA
 CCCGTCGTCATGTTTCGAGGACTCTGTGAATCTGTGGTGAAGCAGACCTCGTGAAGCGCTGAAAAAT
 TTGGGACAATATGCTATGTGATGATGCCATTTAAACGACAGGCTCTAGTGAATTTGAAACATAGA
 TAGTGCCAAAGAATGTGTGACATTTGCTGCAGATGAACCCGTGTACATTGCTGGTCAACAGGCTTTTTTC
 AACTATTCTACAAGCAAAGGATCACTCGGCCAGGAAATACTGATGATCCATCAGGAGGCAACAAAGTTC
 TTCTGCTCTCAATTCAGAATCCGCTTTATCCAATTACAGTGGATGTTTTATATACTGTATGCAACCTGT
 TGGCAAAGTCAACGATTGTTATATTCAAGAGAAATGGGATACAAGCAATGGTTGAGTTTGAATCAGTC
 CTTTGTGCCAGAAAGCTAAAGCAGCACTCAATGGAGCTGATATATATGCTGGATGTTGCACACTAAAAA
 TTGAATATGCACGGCCAACTCGTCTAAATGTTATTAGGAATGACAATGACAGTTGGGACTACACTAAACC
 ATATTTGGGAAGACGAGATAGAGGAAAGGGTCGCCAGAGACAAGCCATTTGGGAGAACACCTTCTTCG
 TTTAGACATGATGGCTATGGATCCCATGGTCCATTATTGCTTTACCAAGTCGTTACAGAATGGGCTCTC
 GAGATACACCTGAACTTGTGCTTATCCATTACCACAGGCTTCTTCTCTTACATGCATGGAGGAAATCC
 CTCTGGTTCAAGTTGAATGGTTAGTGGATTACATCAACTAAAAATGAATTGTTCAAGAGTCTTCAACCTG
 TTCTGCTTATATGGAATATTGAGAAGGTAATTTATGAAGACCATTCTGGTACAGCACTGGTAGAAA
 TGGGTGATGAGTATGCTGTAGAAAGAGCTGTACACACCTTAATAATGTCAAATTTTGGGAAAAGACT
 TAATGTTTGCCTGTCTAAACAACATTCAGTTGTTCCAAGTCAAATATTTGAGCTGGAGGATGGTACCAGC
 AGCTACAAAGATTTTGAATGAGCAAAAAATAATCGCTTTACAAGTCTGGCCAAGCATCTAAGAATATAA
 TCCAGCCACCCTCTGTGTTTTGCATTATTATAATGTTCCATTGTGTGTACAGAAGAGACCTTCACAAA
 GTTGTGTAATGACCATGAAGTTCTTACATTCATCAATATAAGTGTGTTGATGCAAAACCTTCAGCCAAA
 ACACCTTTCTGGGCTATTAGAATGGGAGTGCAAACTGATGCAGTAGAAGCCCTTACGGCACTGAATCACT
 ATCAGATAAGAGTGCCGAATGGTTCCAATCCCTATACATTGAAGCTTTGCTTTCTACATCATCCCATTT
 A

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC226699 representing NM_001142650
 Red=Cloning site Green=Tags(s)


MLARETYEEDREYESQAKRLKTEEGEIDYSAEEGENRREATPRGGDGGGGGRSFSQPEAGGSHHKVSVS
 PVVHVRGLCESVVEADLVEALEKFGTICYVMMMPFKRQALVEFENIDSAKECVTFAADEPVYIAGQQAFF
 NYSTSKRITRPGNTDDPSGGNKVLLLSIQNPLYPTITVDVLYTVCNPVGKVQRIVIFKRNGIQAMVEFESV
 LCAQKAKAALNGADIYAGCCTLKIEYARPTRLNVIRNDNDSWDYTKPYLGRRDRGKGRQRQAILGEHPSS
 FRHDGYGSHGPLPLPSRYRMGSRDTPELVAYPLPQASSSYMHGNGPSGVVMVSGLHQLKMNCRSRVFNL
 FCLYGNIEKVKFMKTIPGTALVEMGDEYAVERAETHLNNVKLFGKRLNVCVSKQHSVVPSSQIFELEDGTS
 SYKDFAMSKNNRFTSAGQASKNIIQPPSCVLHYYNVPLCVTEETF TKLCDNHEVLTFIKYKVFDAKPSAK
 TLSGLLEWECKTDAVEALTALNHYQIRVPNGSNPYTLKLCFSTSSHL

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8066_b04.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001142650

ORF Size: 1611 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001142650.1](#), [NP_001136122.1](#)

RefSeq ORF: 1614 bp

Locus ID: 92906


UniProt ID: [Q8WVV9](#)

Cytogenetics: 2p22.1

MW: 59.5 kDa

Gene Summary: HNRNPLL is a master regulator of activation-induced alternative splicing in T cells. In particular, it alters splicing of CD45 (PTPRC; MIM 151460), a tyrosine phosphatase essential for T-cell development and activation (Oberdoerffer et al., 2008 [PubMed 18669861]).[supplied by OMIM, Aug 2008]

Product images:



Circular map for RC226699

Western validation with an anti-DDK antibody
 * L: Control HEK293 lysate R: Over-expression lysate