

Product datasheet for **RC226701**

ZFP92 (NM_001136273) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	ZFP92
Synonyms:	ZNF897
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC226701 representing NM_001136273
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAGCCATTCTCCTGACCACGAGACCCAAGGTGCCAGTATCTTTGAGGATGTGTCCGTGTA
CTTCAAAAGACAGAAATGGAAGCTTCTGGACCTCAGACAAAAGGTCTCTACAAGCGGGTGTGCTGGAGAACTA
TAGCCATTTGGTGTCACTGGGATTTTCTTCTCCAAGCCTCACCTGATCTCCCAATTGGAACGAGGGGAA
GGACCTTGGGTAGCAGACATCCCCAGAACCTGGGCCACCGCAGGATTGCACATAGGTGACAGAACACAGA
GCAAGACGTCGACTTCAACGCAGAAAGCATTCTGGACGACAACTCCCCGGGGCCGATCCACAAGGTGGCAA
GGAGGGGCAGGCGCGAGGTCGTCTGTGCTCCAGAGAGGTGCCAGGGCTTGGGGCAGAGTTTCGGCTGCG
GGGCCGAGGGCCCAAGGCGCGGAGAAGCGGTACCTGTGCCAGCAGTGTGGGAAGGCCTTCAGCCGCA
GCTCCAACCTCATCAAGCACCGCATCATCCACAGTGGCGAGAAGCCTTACGCGTGCCCCGAGTGCGGCAA
GCTGTTTCGCCGAGCTTCGCGCTCCTGGAGCACCAGCGCATCCACAGCGGCGAGAAGCCCTACGCCTGC
CCCCGAGTGACAGAACCTTCACGCGCAGCTCCAACTCATCAAGCACCAGGTATCCACAGCGGCGAGC
GGCCCTTCGCTGCGGCGACTGCGGCAACTGTCCGCGCAGCTTCGCGCTCCTGGAGCACGCGCGCT
GCACAGCGGCGAGCGGCCCTACGCGTGCCAGAGTGGCGCAAGGCCTTCAGCCGAGCTCCAACCTCATC
GAGCACAGCGCACGACCGCGCGGAGAGCCCTACGCTGCGGCCAGTGCGCCAAGGCCTTCAAGGGCG
TCTCGAGCTCATCCACCACGCGCAGCCACAGCGGCGAGCGGCCCTTCGCGTGCCGCGAGTGCGGCAA
GGCCTTCGCTGGCGTTTCGGGCTCAGCCAGCACCGGCGCGTGACAGCGGTGAGAAGCCCTACGAGTGC
AGCGACTGCGGCAAGGCCTTCGGCCGGCGGCCAACCTATTCAAGCACCAGGAGTGACAGCGCGCCAGGC
GCCCTGCGAAGGCGGAGACGGCGCGGAGGCTAGCGGGCCCTGGGAGCACCAGGCCCTGGGAGCGCGGTGGC
GGCCACCAGCCCCCGCGGCCGAGCACAGCCGCGAGGCCTTCAGGCCAGCCGCCG

AG**CGGACCG**ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA




Protein Sequence: >RC226701 representing NM_001136273
Red=Cloning site Green=Tags(s)

MAAILLTTRPKVPVSFEDVSVYFTKTEWKLLDLRQKVLKRVMLNYSHLVSLGFSFSKPHLISQLERGE
GPWVADIPRTWATAGLHIGDRTQSKTSTSTQKHSGRQLPGADPQGGKEGQAARSSVLQGAQGLGQSSAA
GPQGPKGAEKRYLCQQCGKAFSRSSNLIKHRIIHSGEKPYACPECGKLFRRSFALLEHQRIHSGEKPYAC
PECSKTFTRSSNLIKHQVIHSGERPFACGDCGKLFRRSFALLEHARVHSGERPYACPECGKAFSRSSNL
EHQRTHRGEKPYACGQCAKAFKGVSQLIHHQRSHSGERPFACRECGKAFRGRSGLSQHRRVHSGEKPYEC
SDCGKAFGRANLKFHQAVHGARRPAKAETARRLAGPGSTGPGSAVAATSPPRPSTAARPSRPSRR

SGPTRRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_001136273

ORF Size: 1248 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001136273.1 , NP_001129745.1
RefSeq ORF:	1251 bp
Locus ID:	139735
UniProt ID:	A6NM28
Cytogenetics:	Xq28
MW:	45.6 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC226701