

Product datasheet for **RC226712**

CCDC121 (NM_001142683) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	CCDC121
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-Entry (PSI00001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC226712 representing NM_001142683
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGGGTCCCAGGGCAGGGCTCACACAGCAAGAAAACGAGGAGCATGCCTGTCATTTGAGCCACAGA
GAACGGGGAGCGGAGCCACTGGAGGACCGGCTGCTCGGGCTTATTCGGTAGCCAGGCGGTTAAACAGTT
CAGGGCTGGACAGCCGGGACTGGAGCAGGGTGCAGTCTCCAGGGTTGCTGGGCAGCACCGAGACCCCTT
GAGCACCGAACGAATAAATACTACGGGAGCTTCCACACTTGCACATTGTTCCCGGAGTTGCAGACGCAGG
TTCCTGATGCTAGCGCTCATTCCTTGGCAGTCACCCCTCAGTGAACACAGATTGCCGTGACCTTCAGGA
TGAATGCTTGGATTCCAGGTTTGTGAGGGCCCCAGAAGGCTCCTCCACCGTATCATAGTCTAATAAAT
AATTTTGTCAAGCCAGAGAAGCTAACAAGGTAGAGACAAGGCTTAAAGAAAAGATAGTGGCGGAAATGA
CGGATCTGAACAAGCATATAAAACAAGCTCAAACCCAGCGGAAACAGCTACTGGAGGAATCCAGGGAGCT
ACACCGAGAAAAGTTACTTGTCCAGGCTGAAAACAGATTCTTCTGGAATACCTGACTAACAAAAGTAA
GAGTACACAGAGCAACCTGAGAAGGTATGGAACAGCTATTTACAAAAAGTGGAGAGATTGAACGAAGAA
GACAAGAAATCAGCCTCCAGATATGCAGAACAAATTCAGTGCCTTAAACAGCGCTCTTGCAAAAGGAAAA
TATCCAATCCAGTTTGAAGCGGAAGTTGCAGGCAATGAGGGACATTGCTATATTAAGGAAAAGCAGGAG
AAAGAAATACAGACATTACAGGAGGAGACAAAGAAAGTCCAAGCTGAGACAGCTTCAAAGACACGGGAAG
TACAGGCCAGCTCCTCCAGGAGAAAGATTACTGGAGAAACAAGTGGAGCGAGCCAGACAGGAGGCTACT
GGGAAAGAGAAAAAGAGAGCTTAATATGAAGGCCAGGCTTGAAGTTGGCAGCAAAAGCGGTTTATT
TTTGAATACTCCTGTGGCATCAACAGAGAGAACCAGCAGTTCAAGAAGGAATTACTGCAGCTAATTGAGC
AAGCCAGAACTAACGGCTACTCAAAGCCACTAGAAAACAGGAAGCAGCAGCTGCAGCAGGAACAGTG
GTATCTGGAGTCCTTAATCCAGGCGAGGCAGAGACTGCAAGGAAGTCATAATCAGTGCCTAAATAGACAG
GATGTTCCAAAGACCACACCCAGTCTCCCAAGGCACCAAATCAAGGATTAATCCAAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA



Protein Sequence:

>RC226712 representing NM_001142683

Red=Cloning site Green=Tags(s)

MGSRGQGSKSKTRSMPIVLSQRTGSGATGGPAARAYSVAEAVKQFRAGPAGTGAGCSLQGCWAAPRPF
EHRTNKLRELSTLAHCSRELQTQVPDASAHSALVTLSELHSCRDLDQDECLDSRFAGPRLPPYHSLIN
NFVKPEKLTKVETRLKEKIVAEMTDLNKHIKQAQTQRKQLLEESRELHREKLLVQAENRFFLEYLTNKTE
EYTEQPEKVWNSYLQKSGEIERRRQESASRYAEQISVLKTALLQENIQSSLKRKLQAMRDIAILKEKQE
KEIQTLEETKKVQAETASKTREVQAQLLQEKRLLEKQLSEPDRLLGKRKRRELNMQAQLKLAARFI
FEYSCGINRENQFKELLQLIEQAQKLATQSHLENRKQQLQQEQWYLESLIQARQLQGSNQLNRQ
DVPKTPPSLPQGTCSRINPK

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Chromatograms:

https://cdn.origene.com/chromatograms/mk8060_b05.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_001142683

ORF Size:

1320 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001142683.2](#)

RefSeq ORF: 1323 bp


Locus ID: 79635

UniProt ID: [Q6ZUS5](#)

Cytogenetics: 2p23.3

MW: 50.4 kDa

Product images:



Circular map for RC226712