

Product datasheet for RC226717L3

ATXN7L1 (NM_138495) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ATXN7L1 (NM_138495) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	ATXN7L1
Synonyms:	ATXN7L4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226717).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_138495
ORF Size:	2214 bp



[View online »](#)

