

Product datasheet for **RC226719**

MYADML2 (NM_001145113) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	MYADML2
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-Entry (PSI00001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC226719 representing NM_001145113
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGGCAGCACCATGGAGCCCCCTGGGGGTGCGTACCTGCACCTGGGCGCCGTGACATCCCCTGTGGGCA
CAGCCCGCGTGTGCAGCTGGCCTTTGGCTGCACTACCTTCAGCCTGGTGGCCACCGGGTGGCTTTGC
GGGCGTCCAGGGCACCTTCTGCATGGCCGCCTGGGGCTTCTGCTTCGCCGTCTCTGCGCTGGTGGTGGCC
TGTGAGTTCACACGGCTCCACGGCTGCCTGCGGCTCTCCTGGGGCAACTTCACCGCCGCCTTCGCCATGC
TGGCCACCCTGCTATGCGCGACGGCTGCGGTCTGTATCCGCTGTACTTTGCCCGCGGGAGTGTCCCC
CGAGCCCGCCGGCTGTGCTGCCAGGGAATTCGCGCTGGCAGCCAGTGTCTTCGCCGGGCTCCTTTCCTG
GCCTACGCTGTGGAGGTGGCCCTGACGCGGGCCCGGCCCGGCCAGGTGAGCAGCTATATGGCCACGGTGT
CGGGGCTCCTCAAGATCGTCCAGGCTTCGTGGCTGCATCATCTTCGGGGCGTGGTCCATGACAGCCG
CTACGGGCGCTACGTGGCCACCCAGTGGTGCCTGGCCGTCTACAGCCTGTGCTTCTGGCCACAGTGGCC
GTGGTGGCCCTGAGTGTGATGGGCCACACAGGGGCTGGGCTGCCCTTTGACCGGCTGGTGGTGGTGT
ACACCTTCCTGGCTGTGCTCCTGTACCTCAGCGCCGCGGTGATCTGGCCAGTCTTCTGTTTCGATCCAA
GTACGGTGAGCCCAACGCGCCCCCACTGTGCTCGGGGCGAGCTGTCCCTGGGACAGCCAGCTGGTGGT
GCCATCTTCACCTACGTCAACCTGCTCCTGTACGTCGTTGACCTCGCCTACTCCAGAGGATTGCTTCG
TGCCAGCCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RC226719 representing NM_001145113
Red=Cloning site Green=Tags(s)

MGSTMEPPGGAYLHLGAVTSPVGTARVLQLAFGCTTFLVAHRRGGFAGVQGTFCMAAWGFCFAVSALVVA
CEFTRLHGCLRRLSWGNTAAFAMLATLLCATAAVLYPLFYARRECSPEPAGCAARDFLAASVFAGLLFL
AYAVEVALTRAPGQVSSYMATVSGLLKIQAFVACIIFGALVHDSRYGRYVATQWCVAVYSLCLFATVA
VVALSVMGGTGGGCGPFDRLVVYTFLLAVLLYLSAAVIWPFVCFDPKYGEKRPNPCARGSCPWDSQLVV
AIFTYNNLLLYVLDLAYSORIVFVPSL



TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8117_b04.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001145113

ORF Size: 921 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001145113.3](#)

RefSeq Size: 2469 bp

RefSeq ORF: 924 bp


Locus ID: 255275

UniProt ID: [A6NDP7](#)


Cytogenetics: 17q25.3

MW: 33 kDa

Product images:



Circular map for RC226719



Western blot validation of overexpression lysate (Cat# [LY428699]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC226719 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).