

Product datasheet for **RC226723**

SLCO1C1 (NM_001145944) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	SLCO1C1
Synonyms:	OATP-F; OATP1; OATP1C1; OATP14; OATPF; OATPRP5; SLC21A14
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence: >RC226723 representing NM_001145944
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTGCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGGAGTTGGAACACTGCTCATTGCAATGCCTCAGTTCTTCATGGAGCAGTACAAATATGAGAGATATT
 CTCCTTCCTCCAATTCCACTCTCAGCATCTCTCCGTGTCTCCTAGAGTCAAGCAGTCAATTACCAGTTTC
 AGTTATGGAAAAATCAAAATCCAAATAAGTAACGAATGTGAAGTGAGACTAGCTCTTCCATGTGGATT
 TATGTTTTCTGGGCAATCTTCTCGTGGGAATAGGAGAAATCCCATTGAGCCTTTGGGCATTGCCTACC
 TGGATGATTTTGGCAGTGAAGACAATGCAGCTTTCTATATTGGGTGTGTGCAGACGGTTGCAATTATAGG
 ACCAATCTTTGGTTTCTGTTAGGCTCATTATGTGCCAAACTATATGTTGACATTGGCTTTGTAAACCTA
 GATCACATAACCATTACCCAAAAGATCCCAAGTGGGTAGGAGCCTGGTGGCTTGGCTATCTAATAGCAG
 GAATCATAAGTCTTCTGCAGCTGTGCCTTTCTGGTATTTACCAAAGAGTTTACCAAGATCCCAAAGTAG
 AGAGGATTCTAATTCTTCTCTGAGAAATCCAAGTTTATTATAGATGATCACACAGACTACCAACACCC
 CAGGGAGAAAAATGCAAAATAATGGAAATGGCAAGAGATTTTCTTCCATCACTGAAGAATCTTTTGGAA
 ACCCAGTATACTTCTATATTTATGTACAAGCACTGTTCAAGTTCAATTCTCTGTTCCGCATGGTGACGTA
 CAAACCAAAGTACATTGAGCAGCAGTATGGACAGTCATCCTCCAGGGCCAACTTTGTATCGGGCTCATC
 AACATTCAGCAGTGGCCCTTGGAAATTTCTCTGGGGGATAGTTATGAAAAATTCAGAATCAGTGTGT
 GTGGAGCTGCAAACTCTACTTGGGATCATCTGTCTTTGGTTACCTCCTATTTCTTTCCCTGTTTGCAC
 TGGCTGTGAAAAATCTGATGTGGCAGGACTAACTGTCTCTACCAAGGAACCAACCTGTCTCTTATCAT
 GAACGAGCTCTCTTTTCAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAG
 GTGAAATGGAATCACATATGATCAGCTTGTCTTGTGGTGTCAAACCTCCAACAGGAGTGGAAAAA
 TATTATATTTTACAAGTCACTTGTGTGGGAATGCAGCTTCTAAATCCGGAATTCCTCAGGCATAGTG
 GGAAGATGTCAGAAAGACAATGGATGTCCCAATGTTTCTGTATTTCTTGAATTTCAAGTCAATCATCAT
 CCTATACTTTATCCCTAGTGGCATACCTGGATACATATTACTTCTGAGGTGCATTAAGCCACAGCTTAA
 GTCTTTTGCCTTGGTATCTACACATTAGCAATAAGAGTCTTGCAGGAATCCAGCTCCAGTGTATTTT
 GGAGTTTGTGATGATACTTCAATGCCTCAAAATGGGGATTTAAAGAGTGTGGAAGTAGAGGATCATGCAGAT
 TATATGATTCAAAATGTCTTCAAGTACCAATCAAAAGCATTCCCGCATCTCATTGCTACAGCATCCAGA
 TTTACACAATGCCACTGACACTAACAAGTTTTCATGTCATTTCAGTCTTGTGAAGACATATATCTGGG
 ACTAAGTGTGATACTGGGCACAGTGTCAATTCTCTAAGCATTGCAGTACTTTTCATTTAAAGAAAAAT
 TATGTTTCAAAACACAGAAGTTTATAACCAAGAGAGAAAGAACAATGGTGTCTACAAGATTCCAAAAG
 AAAATTACACTACAAG

ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC226723 representing NM_001145944
 Red=Cloning site Green=Tags(s)


MGVGTLLIAMPQFFMEQYKYERYSPSSNSTLSISPCLLESSSQLPVSVMESKSKISNECEVDTSSTMWI
 YVFLGNLLRGIGETPIQLGIAYLDDFAEDNAAFYIGCVQTVAIIGPIFGFLGSLCAKLYVDIGFVNL
 DHITITPKDPQWVGAWWLGYLIAGIISLLAAMPFWYLPKSLPRSQSREDSNSSEKSKFIIDHDTDYQTP
 QGENAKIMEMARDFLPSLKNLFGNPVYFLYLCTSTVQFNSLFGMVTYKPKYIEQQYQSSSRANFVIGLI
 NIPAVALGIFSGGIVMKKFRISVCGAAKLYLGSSVFGYLLFLSLFALGCENSVDAGLTVSYQGTGPVSYH
 ERALFSDCNSRCKCSETKWPMCGENGITYVSACLACGQTSNRSKNIIIFYNCTCVGIAASKSGNSSGIV
 GRCQKDNCGPQMFLYFLVISVITSYTLGGLIPGYILLRICKPQLKSFALGIYTLAIRVLGAPVYF
 GVLIDTSLCKWGFKRCSRGSCRLYDSNVFRYQIKSIPASHCYSIPDLHNATDNTNKFSCHFACKTYISG
 TNCDTGHSVNSPKHCSTFFHFEKLCFKTQKFYNQERKNNGVYKIPKGLHYK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8069_g01.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001145944

ORF Size: 1836 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001145944.1, NP_001139416.1

RefSeq ORF: 1839 bp

Locus ID: 53919

UniProt ID: [Q9NYB5](#)


Cytogenetics: 12p12.2

Protein Families: Transmembrane

MW: 67.5 kDa

Gene Summary: This gene encodes a member of the organic anion transporter family. The encoded protein is a transmembrane receptor that mediates the sodium-independent uptake of thyroid hormones in brain tissues. This protein has particularly high affinity for the thyroid hormones thyroxine, tri-iodothyronine and reverse tri-iodothyronine. Polymorphisms in the gene encoding this protein may be associated with fatigue and depression in patients suffering from hyperthyroidism. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2009]

Product images:



Circular map for RC226723