

Product datasheet for **RC226730**

FAM47E (NM_001136570) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	FAM47E
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC226730 representing NM_001136570
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGCGGACCGCAGGCGGCGGCTCCGGCCGGGGACGTTGGCCCCGGTGCGCGAGGGCGTGAAGTGCAGGT
CCAGATGTTTCACAAAGCACAAGAACGGGCTGAAGTCCCCACCTCTCTGCACAGCCGCGAGTTGGTATT
TCCAAGAAAGGGGCTGGACGACTTCAGGAAGGGCTGCCCGCCTTGCACTGGTCTGGTGAAGTCCCT
GTGGAGGGCTTTTCCGCCAGATTTATCACAGAGCTCCCCAACTGGCCCCAAAGAAGAGGCAGATCAAGC
TGCTCAAGGAAGCAGACGTGCTTTCCAAGCTCTCGCCAGCCAGCAGGCTCGGAAGGCATTCTGGAGGA
CGTGGAGGCCACCTGACCCACATCCCTTAGCGCTCTACCTGAATCTGGAAGAAGCTATGCCCATAGAG
CTCTTATCAAAGGTGCTGGAAGTGCTTGATCCTGACCGGAAGCTGGAGGACACATGGGCTTATTGTGAGG
ACACCAGGAAAGGAATGAAGGAACCCACGAAGCTTTTAAAAAACATTCTACCCAAGTCTACCTGGGACC
TTCCAAGAAGACGTCTGTGTCAAACGAGGCCAATGGCTTTATGAAGAAAGCCACATAAAATGGATTTG
CTCCATGAAATGGTCTCGTCTGGTCTTCATGAAATGGAATTTCCGACATCGATGAAGAGTTTCATCT
TGAAACAGTTTGACATTGACTATGAGACCAAAACCAAGCCATGATGCGCTCCACACGATGAAGCTAAATCA
GGTTCCTCTGGAGCTAAAGCGTAGTGTGGGGCTCAGTAACTGCAGGAGACAGAGTTCTTCCAGAACTA
GGCTATGAGAGGAACTCCAGAAACCACAGAATCCTTATAAGCCAAAGTGGGTGAAGATGAGGTATGGAG
CATGGTATTTGAACCCCAAGTTGTGGAAGGCAAGAGTAGACGAGCCTCTGGTTGACCTGAGGTCTC
ACATAAGGCTCAAGAGGAGAATTTTAAAAAGGAGCTGCAGGAACAGGAGGAGTTACTTGCAGACCTTCAC
GGAACAGTTGCCTTTAAGGATTTCTTCTAAGCAGGGGCTACAGGACGCCACGTTTCCTTGAGAATATGT
ATATCGGGAAGGAATGTAACGTGCATGTAATAAGACTCCTATAAAACGAACCTCAAGCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

>RC226730 representing NM_001136570

Red=Cloning site Green=Tags(s)

MADRRRLRPGLAPVREGVNCRSRCFTKHKNGLKFPTSLHSRQLVFPKGLDDFRKGCPPTGLVTQVP
VEGFLPQIYHRAPQLAPKKRQIKLLKEADVL SKLSPAQQARKAFLEDVEAHLTPHPLALYLNLEEAMPIE
LLSKVLEVLDPDRKLEDTWAYCQDTRKGMKEPTKLLKKHSTQVYLGPSKKTSVSNAGQWL YEEKPHKMDL
LHENGPRPGLHENGISDIDEEF ILKQFDIDYETKPSHDALHTMKLNQVPLELKRSVGLSKLQETEFFQKL
GYERKLQKPQNPYPKPKVMRYGAWYLNPKLWKKQRVDEPLVDPEVSHKAQEENFKKELQEQEELLADLH
GTVAFKDFILSRGYRTPRFLENMYIGKECKRACNKTPIKRTQA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms:

https://cdn.origene.com/chromatograms/mk8031_e04.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_001136570

ORF Size:

1179 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001136570.3](#)

RefSeq ORF: 1182 bp


Locus ID: 100129583

UniProt ID: [Q6ZV65](#)


Cytogenetics: 4q21.1

MW: 45.5 kDa

Product images:



Circular map for RC226730



Western blot validation of overexpression lysate (Cat# [LY427926]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC226730 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).