

## Product datasheet for **RC226731L4V**

### IFT81 (NM\_001143779) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	IFT81 (NM_001143779) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IFT81
Synonyms:	CDV-1; CDV-1R; CDV1; CDV1R; DV1; SRTD19
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001143779
ORF Size:	2028 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226731).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001143779.1</a>
RefSeq Size:	2742 bp
RefSeq ORF:	2031 bp
Locus ID:	28981
UniProt ID:	<a href="#">Q8WYA0</a>
Cytogenetics:	12q24.11
MW:	80.2 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene, together with IFT74, forms a tubulin-binding module of intraflagellar transport complex B. This module is involved in transport of tubulin within the cilium, and the encoded protein is required for ciliogenesis. Mutations in this gene are a cause of short-rib polydactyly syndromes. [provided by RefSeq, Dec 2016]