

Product datasheet for **RC226736L3V**

PTPH1 (PTPN3) (NM_001145369) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | PTPH1 (PTPN3) (NM_001145369) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | PTPH1 |
| Synonyms: | PTP-H1; PTPH1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001145369 |
| ORF Size: | 2742 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC226736). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001145369.1 , NP_001138841.1 |
| RefSeq Size: | 8783 bp |
| RefSeq ORF: | 2349 bp |
| Locus ID: | 5774 |
| UniProt ID: | P26045 |
| Cytogenetics: | 9q31.3 |
| Protein Families: | Druggable Genome, Phosphatase |
| MW: | 104 kDa |



[View online »](#)

Gene Summary:

The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This protein contains a C-terminal PTP domain and an N-terminal domain homologous to the band 4.1 superfamily of cytoskeletal-associated proteins. P97, a cell cycle regulator involved in a variety of membrane related functions, has been shown to be a substrate of this PTP. This PTP was also found to interact with, and be regulated by adaptor protein 14-3-3 beta. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2009]