

Product datasheet for **RC226739L3V**

OCC1 (C12orf75) (NM_001145199) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	OCC1
Synonyms:	AGD3; OCC-1; OCC1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001145199
ORF Size:	189 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC226739).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001145199.1
RefSeq Size:	1383 bp
RefSeq ORF:	192 bp
Locus ID:	387882
UniProt ID:	Q8TAD7
Cytogenetics:	12q23.3



MW: 6.4 kDa