

Product datasheet for **RC226743**

C20orf173 (NM_001145350) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	C20orf173
Synonyms:	dJ47704.4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC226743 representing NM_001145350
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGAAGCGCTGGCAGATTTTTGTCCTGTGGGTCTTTGGGTGCTCATCCTCTGGCTGATGACCCCTATC
TGGATCTGACACCTGAATCAGACACCCAGGAAAAACGAATGTAAGTGGTACCACAGCATTGCGACTGCC
TTGGTTCAGTTCGGGAAGTGTGGCTGCCCTTCTGAGACCTCAACTGCTCCTCTGCCACACACAGCC
GACGAATGGAAGTGGCTTGATGCGTGCTCCAGGAAGACTATGGGTACCTGATGAGGACAAGAGAGTCTA
TGACCTCTGACACTGTGCTCTGGTGGTTGGGCATGAACCTCAGGGAGCGAGCTTGGGAAATTGTGGAGGAA
GCTGTTAAAGGGATTCCAGGCTCTCGGTGAGCCATTTGATTTCTATTGTGGGACTTGTGTGCTGTTG
GGGCGCCACAGATCCCGCAGGGCTCCAGCCTTGGCAACGACATCGACCAATACCCCGTGGTTTTACGGA
ATGCCAGCGACAGGGCTCCTGGATGCAGCTGGAGATGCTACTGCGGAAGCTCTCTGACCTGGTGTGGAC
TTCAGATGCTCTAAGTATAAGATTTTGAAGATGGTCTGGTCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC226743 representing NM_001145350
Red=Cloning site Green=Tags(s)

MKRWQIFVLVFWVLILWLMTPYLDLTPESAPQEKRMYLVPQHCDPWFSSGKCGCPSETLNCSSCHHTA
DEWNWLDACSRKTMGYLMRTRESMTSDTVLWLMNSGSELGKLWRKLFKGIPLSVSHDFYCGTCVLL
GRPQIPQSSSLGNDIDQYPVFRNASDQGSWMQLEMLLRKLSLWVTS DALSDKILEDGLVP


TRTRPLEQKLISEEDLAANDILDYKDDDDKV



Chromatograms: https://cdn.origene.com/chromatograms/mk8031_e02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM 001145350

ORF Size: 606 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001145350.2
RefSeq ORF:	609 bp
Locus ID:	140873
UniProt ID:	Q96LM9
Cytogenetics:	20q11.22
MW:	23.2 kDa

Product images:



Circular map for RC226743