

Product datasheet for **RC226743L3V**

C20orf173 (NM_001145350) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	C20orf173
Synonyms:	dJ477O4.4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001145350
ORF Size:	606 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC226743).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001145350.1
RefSeq ORF:	609 bp
Locus ID:	140873
UniProt ID:	Q96LM9
Cytogenetics:	20q11.22
MW:	23.2 kDa

