

Product datasheet for RC226751L4V

EFCAB2 (NM_001143943) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	EFCAB2
Synonyms:	CFAP200; DRC8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001143943
ORF Size:	399 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC226751).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001143943.1](#), [NP_001137415.1](#)

RefSeq ORF: 402 bp

Locus ID: 84288

UniProt ID: [Q5VUJ9](#)

Cytogenetics: 1q44

MW: 15.3 kDa



Gene Summary:

The gene encodes a protein that contains two EF-hand calcium-binding domains although its function has yet to be determined. Alternatively spliced transcripts have been observed.
[provided by RefSeq, Mar 2014]