

## Product datasheet for **RC226766A1V**

### Human ATXN7L3B (NM\_001136262) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human ATXN7L3B (NM\_001136262) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** ATXN7L3B  
**Synonyms:** Inc-SCA7  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC226766 representing NM\_001136262  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGAGGAAATTCGTTGGCCAACCTGGATACTAACAAGCTAGAGGCCATCGCTCAGGAGATTTACGTAG  
ACCTGATAGAGGATTCTGTTTGGGATTCTGCTTTGAGGTGCACCGGCCAGTCAAGTGTGGCTACTTCTA  
CCTGGAGTTCGCAGAGACTGGTAGCGTGAAGGATTTGGCATTAGCCAGTGAAGACAAAGGAGCGTGC  
CGCCTCCCGCTTTGCTCCCTTCCCGGAGAACCTGGGAATGGGCCTGATCAGCAGCTGCAGCGCTCACCTC  
CGGAATCCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC226766 representing NM\_001136262  
**Red**=Cloning site **Green**=Tags(s)

MEEISLANLDTNKLEAIAQEIYVDLIEDSCLGFCFEVHRAVKCGYFYLEFAETGSVKDFGIQPVEDKGAC  
RLPLCSLPGEPGNGPDQQLQRSPPPEFQ

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_001136262  
**ORF Size:** 291 bp



[View online »](#)

<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u>NM_001136262.1</u>
<b>RefSeq Size:</b>	3701 bp
<b>RefSeq ORF:</b>	294 bp
<b>Locus ID:</b>	552889
<b>UniProt ID:</b>	<u>Q96GX2</u>
<b>Cytogenetics:</b>	12q21.1
<b>MW:</b>	10.8 kDa