

Product datasheet for **RC226774**

SLCO1C1 (NM_001145945) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	SLCO1C1
Synonyms:	OATP-F; OATP1; OATP1C1; OATP14; OATPF; OATPRP5; SLC21A14
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >RC226774 representing NM_001145945
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGACACTTCATCCAAAGAAATATCCAGTTGTTCTGCAAACTTCAGTGCAACCTGTTGGAAGGCCTT
 CTTTTAAACAGAAATATCCCTCCTCAGAAGAAAAGCAACCATGCTGTGGTGAACATAAGGTGTTCTGTG
 TGCCTTGCTTTTGTCTTACTTTGCCAAGCATTGGCAGAAGGCTATCTGAAGAGCACCATCACTCAGATA
 GAGAGAAGGTTTGATATCCCTTCTTCACTGGTGGGAGTTATTGATGGTAGTTTTGAAATTGGGAATCTCT
 TAGTTATAACATTTGTTAGCTACTTTGGAGCCAACTTCACAGGCCAAAATAATTGGAGCAGGGTGTGT
 AATCATGGGAGTTGGAACACTGCTCATTGCAATGCCTCAGTTCTTTCATGGAGCAGTACAAATATGAGAGA
 TATTCTCCTTCTCCAATTCCACTCTCAGCATCTCTCGTGTCTCCTAGAGTCAAGCAGTCAATTACCAG
 TTTCAAGTTATGAAAAATCAAAATCCAAATAAGTAACGGGTGTGTGCAGACGGTTCGAATTATAGGACC
 AATCTTTGGTTTCTGTAGGCTCATTATGTGCCAACTATATGTTGACATTGGCTTTGTAACCTAGAT
 CACATAACCATACCCCAAGATCCCAAGTGGTAGGAGCCTGGTGGCTTGGCTATCTAATAGCAGGAA
 TCATAAGTCTTCTGACAGTGTGCCTTTCTGGTATTTACCAAGAGTTTACCAAGATCCCAAGTAGAGA
 GGATTTCAATTTCTCTCTGAGAAATCCAAGTTTATTATAGATGATCACACAGACTACCAACACCCAG
 GGAGAAAATGCAAAAATATGGAAATGGCAAGAGATTTTCTCCATCACTGAAGAATCTTTTGGAAACC
 CAGTATACCTTCTATATTTATGTACAAGCACTGTTCAAGTTCAATTCTGTTTCGGCATGGTGACGTACAA
 ACCAAAGTACATTGAGCAGCAGTATGGACAGTCATCCTCCAGGGCCAACTTTGTGATCGGGCTCATCAAC
 ATTCCAGCAGTGGCCTTGAATATTTCTGTTGGGGGATAGTTATGAAAAATTCAGAATCAGTGTGTGTG
 GAGCTGCAAACTCTACTTGGGATCATCTGTCTTTGGTTACCTCTATTTCTTTCCCTGTTTGCAGTGGG
 CTGTGAAAATTTCTGATGTGGCAGGACTAACTGTCTCCTACCAAGGAACCAAACTGTCTCTTATCATGAA
 CGAGCTCTCTTTTCAAGTTCGAATCAAGATGCAATGTTTCAGAGACAAAATGGGAACCCATGTGCGGTG
 AAAATGGAATCACATATGTATCAGCTTGTCTTGTCTGTTGTCAAACCTCCAACAGGAGTGGAAAAATAT
 TATATTTTCAACTGCCTTGTGTGGGAATTCAGCTTCTAAATCCGGAATTCCTCAGGCATAGTGGGA
 AGATGTGAGAAAGACAATGGATGTCCCAATGTTTCTGTATTTCTTGTAAATTCAGTCATCACATCCT
 ATACTTTATCCCTAGGTGGCATACCTGGATACATATTACTTCTGAGGTGCATTAAGCCACAGCTTAAGTC
 TTTTGCCTTGGGTATCTACACATTAGCAATAAGAGTTCTTGCAGGAATCCAGCTCCAGTGTATTTTGA
 GTTTTGATTGATACTTCATGCCTCAATGGGGATTTAAAGATGTGGAAGTAGAGGATCATGCAGATTAT
 ATGATTCAAATGTCTTACAGCATATATATCTGGGACTAACTGTGATACTGGGCACAGTGTCAATTCCT
 AAGCATTGCAGTACTTTTCATTTTAAAGAAAAATTATGTTTCAAAACACAGAAGTTTATAACCAAGAGA
 GAAAGAACATGGTGTCTACAAGATCCAAAGGAAAATTACACTACAAGTGATCATCTGCTACAACCCA
 ACTACTGGCCAGGCAAGGAACTCAACTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC226774 representing NM_001145945
 Red=Cloning site Green=Tags(s)


MDTSSKENIQLFCKTSVQPVGRPSFKTEYPSSSEKQPCCGELKVFLCALSFVYFAKALAEGLKSTITQI
 ERRFDIPSSLVGVIDGSFEIGNLLVITFVSYFGAKLHRPKIIGAGCVIMGVGTLIIAMPQFFMEQKYER
 YSPSSNSTLSISPLLESSQLPVSVMEKSKISNGCVQTVAIIGPIFGFLLGSLCAKLYVDIGFVNLD
 HITITPKDPQWVGAWWLYLIAGIISLLAAMPFWYLPKSLPRSQSREDSNSSSEKSKFIIDDHTDYQTPQ
 GENAKIMEMARDFLPSLKNLFGNPVYFLYLCTSTVQFNSLFGMVTYKPKYIEQQYQSSSRANFVIGLIN
 IPAVALGIFSGGIWMKKFRISVCGAAKLYLGSSVFGYLLFLSLFALGCNSDVAGLTVSYQGTKPVSYHE
 RALFSDCNSRCKCSETKWEPMCGENGITYVSACLQAGCQTSNRSGKNIIFYNCTCVGIAASKSGNSSGIVG
 RCQKDNCGPQMFLYFLVISVITSYLSLGGIPGYILLRLCIKPKLKSFALGIYTLAIRVLGIPAPVYFG
 VLIDTSLKLGWFKRCGSRGSCRLYDSNVFRHIYLGTLVILGTVSILLIAVLFILKKNYVSKHRSFITR
 ERTMVSTRFQKENYTTSDHLLQPNYWPGETQL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8006_g12.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001145945

ORF Size: 1989 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001145945.1](#), [NP_001139417.1](#)

RefSeq ORF: 1992 bp

Locus ID: 53919

UniProt ID: [Q9NYB5](#)


Cytogenetics: 12p12.2

Protein Families: Transmembrane

MW: 73.1 kDa

Gene Summary: This gene encodes a member of the organic anion transporter family. The encoded protein is a transmembrane receptor that mediates the sodium-independent uptake of thyroid hormones in brain tissues. This protein has particularly high affinity for the thyroid hormones thyroxine, tri-iodothyronine and reverse tri-iodothyronine. Polymorphisms in the gene encoding this protein may be associated with fatigue and depression in patients suffering from hyperthyroidism. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2009]

Product images:



Circular map for RC226774