

Product datasheet for **RC226779L3V**

EB2 (MAPRE2) (NM_001143827) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	EB2
Synonyms:	CSCSC2; EBI; EB2; RP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PSI00092)
Tag:	Myc-DDK
ACCN:	NM_001143827
ORF Size:	945 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC226779).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001143827.1](#)

RefSeq ORF: 948 bp

Locus ID: 10982

UniProt ID: [Q15555](#)

Cytogenetics: 18q12.1-q12.2

Protein Families: Druggable Genome



MW: 35.6 kDa

Gene Summary: The protein encoded by this gene shares significant homology to the adenomatous polyposis coli (APC) protein-binding EBI gene family. This protein is a microtubule-associated protein that is necessary for spindle symmetry during mitosis. It is thought to play a role in the tumorigenesis of colorectal cancers and the proliferative control of normal cells. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jan 2012]