

Product datasheet for **RC226787L1V**

NEURL1B (NM_001142651) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	NEURL1B
Synonyms:	hNeur2; neur2; NEURL3; RNF67B
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001142651
ORF Size:	1665 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC226787).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001142651.1
RefSeq ORF:	1668 bp
Locus ID:	54492
UniProt ID:	A8MQ27
Cytogenetics:	5q35.1
MW:	59.1 kDa



Gene Summary:

E3 ubiquitin-protein ligase involved in regulation of the Notch pathway through influencing the stability and activity of several Notch ligands.[UniProtKB/Swiss-Prot Function]