

## Product datasheet for **RC226791**

### FLYWCH2 (NM\_001142500) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** FLYWCH2 (NM\_001142500) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** FLYWCH2  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RC226791 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCCCTGCCGAGCCAGCGAGCAGGAGGGTGAAGGTCAGCCAGGAGCCATCCCCAAGC  
CAGGCACAGAAGTCATCCCGCAGCCCCAGGAAGCCAGAAAGTTCTCCAAACTGGTCTCACAGC  
CTCCAAAGACAGCACCAAGGTGGCGGGGGCCAAGCGCAAGGGTGTGACTGTGTCTATGTCCTGGGGTG  
CCCGCCCCGCCACCTTGCCAAGGCCCTCCTCCAGACCCACCCCGAGGCCAGCGGGCCATTGAGGCAG  
CCCCTCAGGAGCCTGAGCAGAAACGGAGCAGGCCAGGCCAGCACAGACAGAACAAGACAGTGGATT  
AGCAGCGGGCCCTCTGAGGCTGCTGGGAGAAGTTTGCCCCCTGCTCTGTGGCGCCCGCAAGTCCCTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC226791 protein sequence  
Red=Cloning site Green=Tags(s)  
MPLPEPSEQEGESVKASQEPSPKPGTEVIPAAPRKPRKFSKLVLLTASKDSTKVAGAKRKGVHCVMSLGV  
PGPATLAKALLQTHPEAQRRAIEAAPQEPEQKRSRQDPGTDRTEDSGLAAGPPEAAGENFAPCSVAPGKSL

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk6431\\_a07.zip](https://cdn.origene.com/chromatograms/mk6431_a07.zip)

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_001142500

**ORF Size:** 420 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


**RefSeq:** [NM\\_001142500.1](#), [NP\\_001135972.1](#)

**RefSeq Size:** 1040 bp


**RefSeq ORF:** 423 bp

**Locus ID:** 114984


UniProt ID: [Q96CP2](#)  
 Cytogenetics: 16p13.3  
 MW: 14.6 kDa

**Product images:**


Circular map for RC226791



Western blot validation of overexpression lysate (Cat# [LY428134]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC226791 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified FLYWCH2 protein (Cat# [TP326791]). The protein was produced from HEK293T cells transfected with FLYWCH2 cDNA clone (Cat# RC226791) using MegaTran 2.0 (Cat# [TT210002]).