

Product datasheet for **RC226802**

Glutathione S Transferase kappa 1 (GSTK1) (NM_001143680) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Glutathione S Transferase kappa 1 (GSTK1) (NM_001143680) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	GSTK1
Synonyms:	GST; GST13; GST 13-13; GST13-13; GSTK1-1; hGSTK1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC226802 representing NM_001143680 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGGGCCCTGCCGCGCACCGTGGAGCTTTCTATGACGTGCTGTCCCCTACTCCTGGCTGGGCTTCG
AGATCCTGTGCCGTATCAGAATATCTGGAACATCAACCTGCAGTTGCGGCCAGCCTCATAACAGGGAT
CATGAAAGACAGTGGAACAAGCCTCCAGGTCTGCTTCCCGCAAAGGACTATACATGGCAAATGACTTA
AAGCTCCTGAGACACCATCTCCAGATTCATCCACTTCCCAAGGATTTCTTGTCTGTGATGCTTGAAA
AAGGAAGTTTGTCTGCCATGCGTTTCTCACCGCCGTAAGTGGAGCATCCAGAGATGCTGGAGAAAGC
GTCCCGGGAGCTGTGGATGCGCGTCTGGTCAAGGGCTGCAGAGAAGGCTGGTATGTCTGCAGAACAGCC
CAGGGACTTCTGGAAAAGATCGCAACGCCAAAGGTGAAGAACCAGCTCAAGGAGACCACTGAGGCAGCCT
GCAGATACGGAGCCTTTGGGCTGCCATCACCGTGGCCATGTGGATGGCCAAACCCACATGTTATTTGG
CTCTGACCGGATGGAGCTGCTGGCGCACCTGCTGGGAGAGAAGTGGATGGGCCCTATACCTCCAGCCGTG
AATGCCAGACTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC226802 representing NM_001143680
 Red=Cloning site Green=Tags(s)


MGPLPRTVELFYDVLSPYSWLGFEILCRYQNIWNINLQLRPSLITGIMKDSGNKPPGLLPRKGLYMANDL
 KLLRHHLQIPIHFPKDFLSVMLEKGSLSAMRFLTAVNLEHPEMLEKASRELWMRVWSRAAEKAGMSAEQA
 QGLLEKIATPKVKNQLKETTEAACRYGAFGLPITVAHVDGQTHMLFGSDRMELLAHLLGEKWMGPIPPAAV
 NARL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001143680

ORF Size: 642 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001143680.1 , NP_001137152.1
RefSeq ORF:	645 bp
Locus ID:	373156
UniProt ID:	Q9Y2Q3
Cytogenetics:	7q34
Protein Pathways:	Drug metabolism - cytochrome P450, Glutathione metabolism, Metabolism of xenobiotics by cytochrome P450
MW:	24 kDa
Gene Summary:	This gene encodes a member of the kappa class of the glutathione transferase superfamily of enzymes that function in cellular detoxification. The encoded protein is localized to the peroxisome and catalyzes the conjugation of glutathione to a wide range of hydrophobic substates facilitating the removal of these compounds from cells. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Jan 2009]