

Product datasheet for **RC226804L3V**

NCRI (NM_001145457) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	NCRI
Synonyms:	CD335; LY94; NK-p46; NKP46
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001145457
ORF Size:	909 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC226804).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001145457.1](#), [NP_001138929.1](#)

RefSeq ORF: 912 bp

Locus ID: 9437

UniProt ID: [O76036](#)

Cytogenetics: 19q13.42

Protein Families: Druggable Genome, Transmembrane



Protein Pathways: Natural killer cell mediated cytotoxicity

MW: 34.2 kDa

Gene Summary: Cytotoxicity-activating receptor that may contribute to the increased efficiency of activated natural killer (NK) cells to mediate tumor cell lysis.[UniProtKB/Swiss-Prot Function]