

## Product datasheet for **RC226820L1V**

### SLCO2B1 (NM\_001145211) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLCO2B1 (NM_001145211) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLCO2B1
Synonyms:	OATP-B; OATP2B1; OATPB; SLC21A9
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001145211
ORF Size:	2061 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226820).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001145211.3</a>
RefSeq Size:	4156 bp
RefSeq ORF:	2064 bp
Locus ID:	11309
UniProt ID:	<a href="#">O94956</a>
Cytogenetics:	11q13.4
Protein Families:	Transmembrane
MW:	74.4 kDa



[View online »](#)

**Gene Summary:**

This locus encodes a member of the organic anion-transporting polypeptide family of membrane proteins. The protein encoded by this locus may function in regulation of placental uptake of sulfated steroids. Alternatively spliced transcript variants have been described. [provided by RefSeq, Nov 2010]