

## Product datasheet for **RC226823L4V**

### STOML3 (NM\_001144033) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | STOML3 (NM_001144033) Human Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | STOML3   |
| Synonyms:                 | Epb7.2l; SRO   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)  |
| Tag:                      | mGFP   |
| ACCN:                     | NM_001144033   |
| ORF Size:                 | 846 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC226823).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_001144033.1</a> , <a href="#">NP_001137505.1</a>  |
| RefSeq ORF:               | 849 bp   |
| Locus ID:                 | 161003   |
| UniProt ID:               | <a href="#">Q8TAV4</a>   |
| Cytogenetics:             | 13q13.3  |
| Protein Families:         | Transmembrane  |
| MW:                       | 30.9 kDa   |



[View online »](#)

**Gene Summary:**

Required for the function of many mechanoreceptors. Modulate mechanotransduction channels and acid-sensing ion channels (ASIC) proteins. Potentiates PIEZO1 and PIEZO2 function by increasing their sensitivity to mechanical stimulations.[UniProtKB/Swiss-Prot Function]