

Product datasheet for **RC226863A1V**

Human GAS2 (NM_001143830) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human GAS2 (NM_001143830) AAV Particle
Tag:	Myc-DDK
Symbol:	GAS2
Synonyms:	GAS-2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC226863 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTGCACTGCTCTGAGCCAAAGGTACGCAGTGGACCTGGCCTCTCTGATATGCATCAGTATAGCCAAT
GGCTAGCCAGCAGACATGAAGCTAATTTGCTACCAATGAAAGAAGATCTGGCCTTGTGGTTAACCAATCT
ATTAGGGAAGGAGATTACAGCAGAACTTTTATGGAGAAGTTGGACAATGGTGCCTTGCTCTGTCAACTT
GCAGAACTATGCAGGAGAAATTCAAGGAGAGCATGGATGCTAACAAGCCACAAAGAATCTACCGTTGA
AGAAGATCCCATGCAAAACCAAGTGCACCCTCGGGCTCCTTTTTTGGCAGAGACAATACAGCAAAATTTCTT
ATCCTGGTGCCGAGATTTAGGGGTGGATGAAACGTGTCTATTTGAATCGGAAGGTTTGGTCTCCACAAG
CAACCCAGAGAAGTGTGTCTCTGTCTGTAGAGCTTGGCCGGATTGCAGCCAGGTATGGTGTGGAGCCTC
CTGTTTGATAAAGCTGGAAAAAGAGATTGAACAAGAAGAAACACTTTCTGCCCTTCTCCTTCACCTTC
TCCTTCATCAAAGTCTTCTGGAAAAAGAGTACAGGAACTTACTGGATGATGCAGTGAACGAATTTCT
GAAGATCCTCCTTGCAAATGCCCAAACAAGTTCTGTGTGGAGCGGCTCTCCAAGGAAGATACCGAGTGG
GAGAAAAGATCCTCTTCATTAGGATGCTGCACAACAACATGTCATGGTCCGTGTGGGAGGAGGCTGGGA
AACTTTTGCAGGGTATTTGTTGAAACACGACCCTGCCGAATGCTGCAGATCTCCCGTGTGGATGGCAAA
ACATCCCTATCAAAGCAAATCTCAACTCTAAAGGACATGAATCCAGATAACTACTTGGTGGTCTCTG
CCAGTTATAAGGCTAAGAAGGAAATTAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC226863 protein sequence Red=Cloning site Green=Tags(s) MCTALSPKVRSGPGLSDMHQYSQWLASRHEANLLPMKEDLALWLTNLLGKEITAETFMEKLDNGALLCQL AETMQEKFKEKESMDANKPTKNLPLKKIPCKTSAPSGSFFARDNTANFLSWCRDLGVDETCLFESEGLVLHK QPREVCLCLLELGRIAARYGVEPPGLIKLEKEIEQEETLSAPSPSPSPSSKSSGKKSTGNLLDDAVKRIS EDPPCKCPNKFVERLSQGRYRVGEKILFIRMLHNKHVMVRVGGGWETFAGYLLKHDPICRMLQISRVDGK TSPIQSKSPTLKDMPDNYLVVSASYKAKKEIK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001143830
ORF Size:	939 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001143830.1
RefSeq Size:	2240 bp
RefSeq ORF:	942 bp
Locus ID:	2620
UniProt ID:	O43903
Cytogenetics:	11p14.3
MW:	34.9 kDa