

Product datasheet for **RC226879L3V**

HNT (NTM) (NM_001144058) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | HNT (NTM) (NM_001144058) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | HNT |
| Synonyms: | CEPU-1; HNT; IGLON2; NTRI |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001144058 |
| ORF Size: | 1065 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC226879). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001144058.1 , NP_001137530.1 |
| RefSeq ORF: | 1068 bp |
| Locus ID: | 50863 |
| UniProt ID: | Q9P121 |
| Cytogenetics: | 11q25 |
| Protein Families: | Transmembrane |
| MW: | 39 kDa |



[View online »](#)

Gene Summary:

This gene encodes a member of the IgLON (LAMP, OBCAM, Ntm) family of immunoglobulin (Ig) domain-containing glycosylphosphatidylinositol (GPI)-anchored cell adhesion molecules. The encoded protein may promote neurite outgrowth and adhesion via a homophilic mechanism. This gene is closely linked to a related family member, opioid binding protein/cell adhesion molecule-like (OPCML), on chromosome 11. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2009]