

Product datasheet for **RC226881**

DYNC111 (NM_001135556) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	DYNC111
Synonyms:	DNC11; DNC1C1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence: >RC226881 representing NM_001135556
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGTCTGACAAAAGTGACTTAAAGCTGAGCTAGAGCGCAAAAAGCAGCGCTTAGCACAGATAAGAGAAG
 AGAAGAAACGGAAGGAAGAGGAGAGGAAAAAGAGGCTGATATGCAGCAGAAGAAAGAACCCGTTCA
 GGACGACTCTGATCTGGATCGCAACGACGAGAGACAGAGGCTTTGCTGCAAAGCATTGGTATCTCACCG
 GAGCCGCTCTAGTCCCAACCCCTATGTCTCCCTCCTCGAAATCAGTGAGCACTCCCAAGTGAAGCTGGAA
 GCCAAGACTCAGGCGATCTGGGGCCATTAAACAGGACCCTGCAGTGGGACACAGACCCCTCAGTGCTCCA
 GCTGCAGTCAGACTCAGAAGTGGAAAGAACTGCATAAACTGGGCGTGTCAAAGGTCAACCAAGTGGAT
 TTCCTGCCAAGGGAAGTAGTGTCTACTCAAAGGAGACCCAGACTCCTCTTGCCACGCATCAGTCTGAAG
 AGGATGAGGAAGATGAGGAAATGGTGAATCTAAAGTTGGCCAGGACTCAGAAGTGGAAATCAGGACAA
 AAAACAGGAAGTGAAGGAAGCCCTCCAAGAGAGTTGACAGAGGAAGAAAACAGCAGATCATTATTCA
 GAGGAATTTCTCATCTTTTTGACCGGACAATACGGGTAAATTGAAAGAGCCCTGGCTGAAAGATTCCGACA
 TCTTTTTGACTACAGCGGCCGAGAGTTAGAGGAAAAAGATGGGGATGTTCAAGCTGGAGCCAATCTTTC
 TTTCAATCGTCAGTTCTATGATGAACATTGTTCCAAGCATCGAGTGGTCACTTGTATGGACTGGTCCCTC
 CAGTACCCTGAGCTGATGGTGGCTTCTTACAACAACATGAAGATGCTCCCCATGAACCAGATGGAGTGG
 CCTTGGTTTGAACATGAAGTTTAAAGAAACCAACCAAGAAATACGTCTTCCACTGTCACTCCTCTGTGAT
 GTCGGTCTGCTTCGCCCCGTTTCCATCCTAACTTGGTGGTTGGTGGGACTTACTCGGGCCAGATTGTCCTC
 TGGGACAATCGCAGTCATCGAAGGACTCCAGTGCAGCGGACACCCTTATCAGCTGCTGCACACACGCATC
 CCGTGTACTGTGTAATGTTGTTGGGACCCAGAATGCTCATAACCTCATCACTGTCTCCACTGATGGCAA
 AATGTGTTCTGGAGCCTGGACATGCTCTCAACTCCACAGGAGAGCATGGAGCTGGTGTACAATAAGTCC
 AAGCCTGTGCTGTTACCGGAATGGCTTCCCAACGGGAGAGCGTCAATAACTTCGTGGTTGGCAGTGAGG
 AAGGTACAGTCTACACGGCTTGTGTCATGGAAGCAAAGCAGGTATTGGTGAGGTCTTTGAAGGTACCA
 AGGGCCAGTGACAGGAATTAAGTCCACATGGCAGTGGGACCAATCGACTTTTCTCACCTGTTTGTGACA
 TCATCATTTGACTGGACTGTCAAAGTGTGGACCAACCAAGCACAACAGCCGCTTACTCCTTTGAAGACA
 ATGCAGACTATGTGTACGATGTGATGGTCCCCCGTGCATCCTGCGCTTTTGCCTGCGTGGACGGGAT
 GGGGCGCTTGGACCTCTGGAACCTCAACAATGACACCGAGGTTCCAACAGCAAGTGTGGCCATTGAGGGG
 GCATCCGCCCTAAACCGTGTCTGTTGGGCCCAAGCTGGCAAAGAAGTTGCTGTTGGGACTCGGAAGGCC
 GTATTTGGGTCTATGACGTTGGAGAGCTTGCAAGTCCCCACAATGATGAATGGACCCGATTGTCAGGAC
 CCTTGTGGAATTCGTGCTAACAGAGCTGATAGCGAGGAGGAAGGCACTGTTGAGTTATCTGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC226881 representing NM_001135556
 Red=Cloning site Green=Tags(s)

MSDKSDLKAELEKQRLAQIREKKRKEEERKKKEADMQQKKEPVQDDSDLDKRRETEALLQSIGISP
 EPPLVPTPMSPKSVSTPSEAGSQDSDLGPLTRTLQWTDPSVLQLQSDSELGRRHLKLGVSQVQVD
 FLPREVVSYSKETQTPLATHQSEDEEDEMVEKVGQDSELENQDKKQEVKEAPPRELTEEEKQIIHS
 EEFLIFFDRTIRVIERALAEDSDIFFDYSGRELEEKDGDVQAGANLSFNRQFYDEHCSKHRVVTMDWSL
 QYPELMVASYNNEADAPHEPDGVALVWNMFKKTTPYVFHCQSSVMSVCFARFHPNLVVGTYSGQIVL
 WDNRSRRTPVQRTPLSAAATHPVYCVNVVGTQNAHNLITVSTDGKMCSWSLDMLSTPQESMELVYNKS
 KPVAVTGMAFPTGDVNNFVVGSEEGTVYACRHGSKAGIGEVFEHQGPVTGINCHMAVGPIDFSHLFVT
 SSFDWTVKLWTTKHNKPLYSFEDNADYVYDVMWSPVHPALFACVDGMGRDLWNLNNDTEVPTASVAIEG
 ASALNRVWAQAGKEVAVGDSEGRIVYVDVGELAVPHNDEWTRFARTLVEIRANRADSEEGTVELSA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001135556

ORF Size: 1884 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001135556.1](#), [NP_001129028.1](#)

RefSeq ORF: 1887 bp

Locus ID: 1780


UniProt ID: [Q14576](#)

Cytogenetics: 7q21.3

MW: 70.6 kDa

Gene Summary: Acts as one of several non-catalytic accessory components of the cytoplasmic dynein 1 complex that are thought to be involved in linking dynein to cargos and to adapter proteins that regulate dynein function. Cytoplasmic dynein 1 acts as a motor for the intracellular retrograde motility of vesicles and organelles along microtubules. The intermediate chains mediate the binding of dynein to dynactin via its 150 kDa component (p150-glued) DCNT1. May play a role in mediating the interaction of cytoplasmic dynein with membranous organelles and kinetochores.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC226881