

Product datasheet for **RC226897**

KIAA0652 (ATG13) (NM_001142673) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	KIAA0652
Synonyms:	KIAA0652; PARATARG8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >RC226897 representing NM_001142673
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGAACTGATCTCAATTCCCAGGACAGAAAGGACCTGGACAAGTTTATTAATTTTTTGCCTCAAGA
 CTGTCCAAGTGATTGTCCAGGCTCGGCTTGGTGAAGATTGCACTCGTTCATCATCTTCTCCAACGGG
 TTCAGATTGGTTCAACTTAGCAATCAAGACATCCCAGAGGTTACACATGAAGCAAAGAAGGCACTGGCA
 GGACAGCTGCCTGCAGTCGGGAGGTCCATGTGTGTGGAGATTTCACTTAAGACTTCTGAGGGAGATTCCA
 TGGAGCTGGAATATGGTGTCTTGAATGAATGAAAGTGTGATAAAGAAATCAAAGTTTCTACACGGT
 GTACAACAGACTGTCAATTGCTGCTGAAGTCCCTTCTTGCTATAAAGTGGTGACACCAAGCTATAGGCTC
 TCCAGGAAACAAGGGCATGAATATGTCATATTACAGGATATATTTGGAGAAGTTCAGCTGAGTGGCT
 TAGGAGAAGGCTTCCAGACAGTTCGTGTTGGGACAGTGGCACCCCTGTGGCACCATCACTCTTCTTG
 TGCTTACAGAAATTAAGTGGCATTGCTGTACAGGCAATTTGAGAGGACCCACCTATCATGGGGATT
 ATTATTGATCACTTGTGGACCGTCCCTATCCCAGCTCCTCTCCCATGCACCCCTGCAATTACAGAACTG
 CTGGTGAGGACACTGGAGTAATATACCCGCTGTAGAAGACTCTCAAGAAGTGTGTACCACTCTTTTTC
 CACCTCCCCACCATCCCAGCTCTCAAGCTCTCGCTTTCCTATCAGCCTGCTGCCCTGGGCGTTGGATCA
 GCTGACCTGGCTTATCCAGTAGTGTGCTGCTGGCTTAAATGCTACACACCCCTACCAGCTGATGGTTC
 CTGGGAAGGAAGGTGGGTACCCCTTGCTCCCAACCAAGCTGTCATGGTACCCAGGCTGACCAAGGAG
 ACTGGCAACCTGCACCCCTTCTGACAGAACCCTGTGCTGCCACACCCCTCCAGTAGTGAGGATACTGAA
 ACCGTATCAACAGCAGTGAGGGACGGGCTCCCTCACGATGCTTGGAGACCATCTTGTCCGAAAAG
 TGGGGGCTTTGTCAACAAACCCATTAACAGGTGACCCTGACGAGTTTGGATATACCCCTTGCCATGTT
 TGCTCCCAAGAATTTGGAGCTGGAGGATACCGATCCAATGGTGAATCCTCCAGATTCCCAGAGACTGAA
 TCTCTCTCCAGGGCAGCTGCACTCAGATGGCTCCAGCGGGGAGCAGTGGCAATACCCATGATGACT
 TTGTTATGATAGACTTTAAACAGCTTTTCTAAGATGACATTCTTCCGATGGACCTGGGGACCTTCTA
 TCGGGAGTTTCAGAACCCACCTCAGCTGAGCAGCTCTCCATAGATATTGGAGCACAGTCCATGGCTGAA
 GACTTGGACTCATTACCAGAGAAGCTGGCTGTGCATGAGAAGAATGCCGCGAGTTTGATGCCTTTGTGG
 AAACCTGCGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC226897 representing NM_001142673
 Red=Cloning site Green=Tags(s)


METDLNSQDRKDLDFIKFFALKTVQVIVQARLGEKICTRSSSSPTGSDWFNLAIKDIPEVTHEAKKALA
 GQLPAVGRSMCVEISLKTSEGDSMELEIWCLEMNEKCDKEIKVSYTVYNRLSLLKSLLAITRVTPAYRL
 SRKQGHEHYILYRIYFGEVQLSGLGEGFQTVRVGTGTPVGTITLSCAYRINLAFMSTRQFERTPPIMGI
 IIDHFVDRPYPPSSPMHPCNYRTAGEDTGVYPSVEDSQEVCTTSFSTSPPSQLSSRLSYQPAALGVGS
 ADLAYPVVFAAGLNATHPHQLMVPKKEGGVPLAPNPVHGTQADQERLATCTPSDRTHCAATPSSSEDTE
 TVSNSSEGRASPHDVLETIFVRKVGAFVNKPINQVTLTSLDIPFAMFAPKNLEEDTDPMVNPPDSPETE
 SPLQGSLSHDSGSSGSGNTHDDFVMIDFKPAFSKDDILPMDLGTFFREFQNPQLSSLIDIGAQSMAE
 DLDSLPEKLAVHEKNVREFDAFVETLQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mg4744_d04.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001142673

ORF Size: 1551 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.


The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001142673.2, NP_001136145.1</u>
RefSeq ORF:	1554 bp
Locus ID:	9776
UniProt ID:	<u>O75143</u>
Cytogenetics:	11p11.2
MW:	56.4 kDa
Gene Summary:	The protein encoded by this gene is an autophagy factor and a target of the TOR kinase signaling pathway. The encoded protein is essential for autophagosome formation and mitophagy. [provided by RefSeq, Oct 2016]

Product images:



Circular map for RC226897



Western blot validation of overexpression lysate (Cat# [LY428237]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC226897 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).