

Product datasheet for RC226900L3

BHLHB9 (NM_001142528) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BHLHB9 (NM_001142528) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	BHLHB9
Synonyms:	GASP3; p60TRP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226900).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_001142528
ORF Size:	1641 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001142528.1 , NP_001136000.1
RefSeq ORF:	1644 bp
Locus ID:	80823
UniProt ID:	Q6PI77
Cytogenetics:	Xq22.1
MW:	60.1 kDa
Gene Summary:	This gene is a member of a gene family which encodes proteins with a basic helix-loop-helix domain. Other members of this gene family encode proteins which function as transcription factors, either enhancing or inhibiting transcription depending on the activity of other DNA binding proteins. The coding region of this gene is located entirely within the terminal exon. The encoded protein may be involved in the survival of neurons (PMID: 15034937). Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Sep 2011]

Product images:



Circular map for RC226900L3