

Product datasheet for **RC226917L3V**

TSEN2 (NM_001145392) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	TSEN2
Synonyms:	PCH2B; SEN2; SEN2L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001145392
ORF Size:	1395 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC226917).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001145392.1 , NP_001138864.1
RefSeq Size:	2150 bp
RefSeq ORF:	1398 bp
Locus ID:	80746
UniProt ID:	Q8NCE0
Cytogenetics:	3p25.2



MW: 53.2 kDa

Gene Summary: This gene encodes one of the subunits of the tRNA splicing endonuclease. This endonuclease catalyzes the first step in RNA splicing which is the removal of introns. Mutations in this gene have been associated with pontocerebellar hypoplasia type 2. A pseudogene has been identified on chromosome 4. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Feb 2009]