

Product datasheet for **RC226919L3V**

ELAVL4 (NM_001144774) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ELAVL4 (NM_001144774) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ELAVL4
Synonyms:	HUD; PNEM
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001144774
ORF Size:	1098 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226919).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001144774.1
RefSeq ORF:	1101 bp
Locus ID:	1996
UniProt ID:	P26378
Cytogenetics:	1p33-p32.3
Protein Families:	Druggable Genome
MW:	40.2 kDa



[View online »](#)

Gene Summary:

May play a role in neuron-specific RNA processing. Protects CDKN1A mRNA from decay by binding to its 3' UTR (By similarity). Binds to AU-rich sequences (AREs) of target mRNAs, including VEGF and FOS mRNA.[UniProtKB/Swiss-Prot Function]