

Product datasheet for **RC226939L3V**

DTHD1 (NM_001136536) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Symbol: DTHD1

Mammalian Cell Selection: Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Tag: Myc-DDK

ACCN: NM_001136536

ORF Size: 1848 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC226939).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001136536.2](#)

RefSeq ORF: 1851 bp

Locus ID: 401124

UniProt ID: [Q6ZMT9](#)

Cytogenetics: 4p14

MW: 69.7 kDa



Gene Summary:

This gene encodes a protein which contains a death domain. Death domain-containing proteins function in signaling pathways and formation of signaling complexes, as well as the apoptosis pathway. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2012]