

## Product datasheet for **RC226940L3V**

### SLC1A5 (NM\_001145144) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC1A5 (NM_001145144) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC1A5
Synonyms:	AAAT; ASCT2; ATBO; M7V1; M7VS1; R16; RDRC
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001145144
ORF Size:	939 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC226940).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001145144.1</a>
RefSeq Size:	1737 bp
RefSeq ORF:	942 bp
Locus ID:	6510
UniProt ID:	<a href="#">Q15758</a>
Cytogenetics:	19q13.32
Protein Families:	Druggable Genome, Transmembrane
MW:	33.7 kDa



[View online »](#)

**Gene Summary:**

The SLC1A5 gene encodes a sodium-dependent neutral amino acid transporter that can act as a receptor for RD114/type D retrovirus (Larriba et al., 2001 [PubMed 11781704]).[supplied by OMIM, Jan 2011]