

Product datasheet for **RC226957L4V**

DDX4 (NM_001136034) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | DDX4 (NM_001136034) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | DDX4 |
| Synonyms: | DEAD (Asp-Glu-Ala-Asp) box polypeptide 4; DEAD/H (Asp-Glu-Ala-Asp/His) box polypeptide 4; MGC111074; OTTHUMP00000122546; VASA |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001136034 |
| ORF Size: | 2172 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC226957). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001136034.1 , NP_001129506.1 |
| RefSeq Size: | 2880 bp |
| RefSeq ORF: | 2174 bp |
| Locus ID: | 54514 |
| Cytogenetics: | 5q11.2 |
| MW: | 79.3 kDa |



[View online »](#)

Gene Summary:

DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. This gene encodes a DEAD box protein, which is a homolog of VASA proteins in *Drosophila* and several other species. The gene is specifically expressed in the germ cell lineage in both sexes and functions in germ cell development. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2009]