

Product datasheet for **RC226971**

C11orf86 (NM_001136485) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C11orf86 (NM_001136485) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: C11orf86
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC226971 representing NM_001136485
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGAACAGGGCTGCGAAGTCAGTCCTTGCAGAGCCCCGACCCTCTATGGAAAGCTGCAGGAGCCCT
GGGGGAGGCCCCAGGAGGGCCAACCTCCGCAGGGCGCTAAGCCTCAGACAGGGGCAAGAGAAGTCCAGGTC
CCAGGGCCTCGAGAGAGGCACAGAAGGGCCAGATGCCACTGCCAGGAGCGGGTGCCGGGGAGCCTGGGG
GACACAGAGCAGCTGATCCAAGCCAGCGAAGAGGCAGCCGGTGGTGGCTGAGGCGGTACCAACAGGTAA
GAAGAAGTGGGAGAGCTTTGTGCCATCTTCCCAGCGTGACTCTGAGTCAGCCGGCCTCCCCG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC226971 representing NM_001136485
Red=Cloning site **Green**=Tags(s)

MGTGLRSQSLREPRPSYGLQEPWGRPQEGQLRRALSLRQGQEKSRSQGLERGTEGPDATAQERVPGSLG
DTEQLIQARRRSRWLRRYQQVRRRWESFVAIFPSVTLSPASP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8032_f10.zip

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001136485

ORF Size: 345 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_001136485.3](#)

RefSeq ORF: 348 bp


Locus ID: 254439

UniProt ID: [A6NJI1](#)
 Cytogenetics: 11q13.2
 MW: 13 kDa


Product images:



Circular map for RC226971



Western blot validation of overexpression lysate (Cat# [LY427891]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC226971 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified C11orf86 protein (Cat# [TP326971]). The protein was produced from HEK293T cells transfected with C11orf86 cDNA clone (Cat# RC226971) using MegaTran 2.0 (Cat# [TT210002]).