

Product datasheet for **RC226977L3V**

C17orf49 (NM_001142799) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	C17orf49
Synonyms:	BAP18; HEPIS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001142799
ORF Size:	414 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC226977).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001142799.1
RefSeq ORF:	417 bp
Locus ID:	124944
UniProt ID:	Q8IXM2
Cytogenetics:	17p13.1
MW:	13.8 kDa



Gene Summary:

Component of chromatin complexes such as the MLL1/MLL and NURF complexes.
[UniProtKB/Swiss-Prot Function]