

Product datasheet for **RC226989L3V**

ZBTB43 (NM_001135776) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	ZBTB43
Synonyms:	ZBTB22B; ZNF-X; ZNF297B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001135776
ORF Size:	1401 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC226989).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001135776.1 , NP_001129248.1
RefSeq Size:	5851 bp
RefSeq ORF:	1404 bp
Locus ID:	23099
UniProt ID:	O43298
Cytogenetics:	9q33.3



Protein Families: Transcription Factors

MW: 52.6 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]