

Product datasheet for RC227006

CMKLR1 (NM_001142343) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CMKLR1 (NM_001142343) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CMKLR1
Synonyms:	CHEMERINR; ChemR23; DEZ; RVER1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC227006 representing NM_001142343 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGAATGGAGGATGAAGATTACAACACTTCCATCAGTTACGGTGATGAATACCCTGATTATTTAGACT
CCATTGTGGTTTTGGAGGACTTATCCCCCTTGAAGCCAGGGTGACCAGGATCTTCTGGTGGTGGTCTA
CAGCATCGTCTGCTTCTCGGGATTCTGGGCAATGGTCTGGTGATCATCATTGCCACCTCAAGATGAAG
AAGACAGTGAACATGGTCTGGTTCCTCAACCTGGCAGTGGCAGATTTCTGTTCAACGTCTTCTCCCAA
TCCATATCACCTATGCCGCCATGGACTACCACTGGGTTTTCGGGACAGCCATGTGCAAGATCAGCAACTT
CCTTCTCATCCACAACATGTTACCAGCGTCTTCTGCTGACCATCATCAGCTCTGACCGCTGCATCTCT
GTGCTCCTCCCTGTCTGGTCCCAGAACCACCGCAGCGTTCGCCTGGCTTACATGGCCTGCATGGTATCT
GGGTCTGGCTTTCTTCTTGAGTTCCTCATCTCTCGTCTTCCGGGACACAGCCAACCTGCATGGGAAAAT
ATCCTGCTTCAACAACCTCAGCCTGTCCACACCTGGGTCTTCTCGTGGCCCACTCACTCCCAAATGGAC
CCTGTGGGGTATAGCCGGCAGATGGTGGTACTGTACCCCGTTCCTCTGTGGCTTCTGGTCCCAGTCC
TCATCATCACAGCTTGTACCTACCATCGTGTGCAAACCTGCAGCGCAACCGCTGGCCAAGACCAAGAA
GCCCTTCAAGATTATTGTGACCATCATATTACCTTCTTCTCTGCTGGTGGCCCTACCACACTCAAC
CTCCTAGAGCTCCACCACACTGCCATGCCTGGCTGTCTTTCAGCCTGGGTTTGGCCCTGGCCACTGCC
TTGCCATTGCCAACAGCTGCATGAACCCATTCTGTATGTTTTTCATGGGTCAGGACTTCAAGAAGTTCAA
GGTGGCCCTCTTCTCTCGCTGGTCAATGCTCTAAGTGAAGATACAGGCCACTTCTCCTACCCAGCCAT
AGAAGCTTACCAAGATGTCATCAATGAATGAGAGGACTTCTATGAATGAGAGGGAGACCGCATGCTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA



[View online »](#)

Protein Sequence: >RC227006 representing NM_001142343
 Red=Cloning site Green=Tags(s)

MRMEDEDYNTSISYGDEYPDYLDLSIVVLEDLSPLEARVTRIFLVVYVIVCFGLGILGNLVIIIATFKMK
 KTVNMVWFLNLAVADFLFNVFLPIHITYAAMDYHWVFGTAMCKISNLLIHNMFTSVFLLTIISSDRICIS
 VLLPVWSQNHRSVRLAYMACMVIWVLAFFLSSPSLVFRDTANLHGKISCFNNFSLSTPGSSSWPHTSQMD
 PVGYSRHMVVTVTRFLCGFLVPLIITACYLTIIVCKLQRNRLAKTKKPKIIVTIIITFFLCWCPYHTLN
 LLELHHTAMPGSVFSLGLPLATALAIANSCMNPILYVFMGQDFKKFKVALFSRLVNALSEDTGHSSYPHS
 RSFTKMSSMNERETSMNERETGML

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001142343

ORF Size: 1119 bp


OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001142343.2</u>
RefSeq ORF:	1122 bp
Locus ID:	1240
UniProt ID:	<u>Q99788</u>
Cytogenetics:	12q23.3
Protein Families:	Druggable Genome, GPCR, Transmembrane
MW:	42.1 kDa
Gene Summary:	<p>Receptor for the chemoattractant adipokine chemerin/RARRES2 and for the omega-3 fatty acid derived molecule resolvin E1. Interaction with RARRES2 induces activation of intracellular signaling molecules, such as SKY, MAPK1/3 (ERK1/2), MAPK14/P38MAPK and PI3K leading to multifunctional effects, like, reduction of immune responses, enhancing of adipogenesis and angiogenesis. Resolvin E1 down-regulates cytokine production in macrophages by reducing the activation of MAPK1/3 (ERK1/2) and NF-kappa-B. Positively regulates adipogenesis and adipocyte metabolism. Acts as a coreceptor for several SIV strains (SIVMAC316, SIVMAC239, SIVMACL7E-FR and SIVSM62A), as well as a primary HIV-1 strain (92UG024-2). [UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for RC227006

