

Product datasheet for RC227014L4

TMEM169 (NM_001142311) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM169 (NM_001142311) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	TMEM169
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227014).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

RefSeq ORF: 894 bp
Locus ID: 92691
UniProt ID: Q96HH4
Cytogenetics: 2q35
Protein Families: Transmembrane
MW: 33.6 kDa