

Product datasheet for **RC227021L4V**

KIAA1432 (RIC1) (NM_001135920) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KIAA1432 (RIC1) (NM_001135920) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KIAA1432
Synonyms:	bA207C16.1; CATIFA; CIP150; KIAA1432
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001135920
ORF Size:	3495 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227021).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001135920.2
RefSeq Size:	3831 bp
RefSeq ORF:	3498 bp
Locus ID:	57589
UniProt ID:	Q4ADV7
Cytogenetics:	9p24.1
MW:	131.2 kDa


[View online »](#)

Gene Summary:

The RIC1-RGP1 complex acts as a guanine nucleotide exchange factor (GEF), which activates RAB6A by exchanging bound GDP for free GTP and may thereby required for efficient fusion of endosome-derived vesicles with the Golgi compartment. The RIC1-RGP1 complex participates in the recycling of mannose-6-phosphate receptors. Required for phosphorylation and localization of GJA1.[UniProtKB/Swiss-Prot Function]