

Product datasheet for **RC227030L3V**

SFRS15 (SCAF4) (NM_001145445) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SFRS15 (SCAF4) (NM_001145445) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SFRS15
Synonyms:	SFRS15; SRA4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001145445
ORF Size:	3375 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC227030).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001145445.1 , NP_001138917.1
RefSeq ORF:	3378 bp
Locus ID:	57466
UniProt ID:	O95104
Cytogenetics:	21q22.11
MW:	123.5 kDa


[View online »](#)

Gene Summary:

This gene likely encodes a member of the arginine/serine-rich splicing factor family. A similar protein in Rat appears to bind the large subunit of RNA polymerase II and provide a link between transcription and pre-mRNA splicing. Alternatively spliced transcript variants have been described. [provided by RefSeq, Feb 2009]